INGEZIM EHRLICHIA VET

Prod Ref: 15.EHR.K8-32

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I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the antigen on a solid support (polystyrene plate). When a sample serum contains specific antibodies against *E.canis*, they will bind to the antigen adsorbed on plate. After washing to eliminate all non-fixed material from the sera

sample, we can detect the presence of dog IgG using a specific peroxidase conjugate.

After addition of the substrate a colorimetric reaction will appear.

In this way the presence of colour means the presence of antibodies against *E.canis* in the dog sera, and the absence of colour the absence of specific antibodies

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **IMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Use a new loop for each sera sample.
8. For each utilization of the kit, control positive and negative serum must be tested in a systematic way.

III. STORAGE OF COMPONENTS

The reagents must be storage between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using a washing solution bottle or an automatic washer or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instruction:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense the washing solution (prepared following instructions), proportionally to the wells used.
- Shake delicately the plate, avoiding the contamination between wells. Turn over the plate brusquely to empty the wells.
- Repeat the process as much times is indicated on the instructions of the kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use.
- Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

The sera samples need to be at **1/100** dilution in sera diluent to be tested. You get it using the calibrated loop and following the instruction

VI. WAY TO USE THE CALIBRATED LOOPS:

- ◆ Use one different calibrated loop for each sera sample.
- ◆ For taking the adequate volume introduce only the final loop into the sample container.
- ◆ The captured volume is deposited by simple introduction of the loop into the appropriate well (previously has been add the sera sample diluent).
- ◆ The dilution is homogenized by soft turns of the loop into the well.

VII. PREPARATION OF REAGENTS:

All the reagents are ready to use, without any previous manipulation. Only the washing solution is 10 x concentrated.

9 parts of distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

- **Washing solution preparation:**
Dilute one part of the concentrate washing solution provided in the kit into

Example: 60 ml. of concentrated washing solution into 540 ml. of distilled or desionized water.

VIII. PROCEDURE:

1. Take out the strips or wells to be used and keep it at least 2 hours at room temperature before starting the test.
2. Add **2 drops** of serum diluent on every well where serum samples are going to be tested. Then and using the calibrated loops, add **1 µl** of serum sample on each well following the previous instructions.
The positive and negative controls are ready to use and is only necessary add **2 drops** of each one on the remainder wells.
Seal the plate and **incubate for 10 minutes at room temperature.**
3. **Wash 5 or more times** following the described procedure.
4. Add **2 drops** of conjugate (RED COLOUR), to each well (controls and serum samples).
Seal the plate **and incubate for 10 minutes at room temperature.**
5. Wash 5 or more times following the described procedure.
6. Add **1 drop** of cromogen solution, to each well and **1 drop** of substrate solution in the same way. Keep the plate for **5 min at room temperature.**
7. Add **2 drops** of stop solution.
8. Read the results.

IX. READING AND RESULT INTERPRETATION

The reading must be done by visual colour turn detection.

Kit validation:

We take the controls like reference.

- **Negative control:** There is no colour turn detection (transparent).
- **Positive control:** There is colour turn detection (blue).

Results Interpretation:

- **Negative samples:** There is no colour turn detection.

- **Positive samples:** There is colour turn detection.

▪ **ATTENTION**

The assay must be used as a **qualitative test** (positive / negative). Nevertheless and because the development of colour is directly proportional to the titre of antibodies, us much intense is the colour, us much titre of antibodies specific of ehrlichia will be.

Developed and manufactured in Spain by:

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INMUNOLOGIA Y GENETICA APLICADA, S.A.