INGEZIM PPA DAS 2.0

Prod Ref: 11.PPA.K2

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I. TECHNICAL BASIS

This kit is based on a double antibody sandwich enzymatic immunoassay (DAS or capture Elisa). We make a brief description of the technique below:

The plate is coated with a specific monoclonal antibody (Mab) to VP72 of the ASFV. In this way when a sample is added into the wells, if they contain antigen, the coated Mabs will catch it. After washing a second Mab to VP72 (specific for a different epitope) conjugated with Peroxidase, is added.

If there is antigen captured into the well, the conjugate (Mab, labelled with peroxidase) will bind it and the peroxidase action will be detected after adding the adequate substrate.

Samples to be used in the assay could be any tissue suspected to be infected. For diagnoses purposes, we recommend to use as sample, blood, spleen or lymphatic nodes (tissues in which viral concentration could be high).

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or kit reagents are being handle.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample. Use only distilled water for the preparation of the reagents.
9. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way
10. Substrate and substrate buffer must be handle with care.

III. STORAGE OF COMPONENTS

All components must be stored between +2°C and +8°C. **Once opened**, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

WARNING! Stop Solution may precipitate, if it happens, warm it at +37°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use.
- Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the 25x concentrate washing solution provided in the kit with 24 parts of distilled or deionized water. When ready this solution remains stable at +4°

- **Controls:**

Before use a 1/50 dilution must be done in diluent. Once opened, controls are stable for 1 month between +2°C and +8°C.

For long term storage, we recommend to make aliquots of the controls and store it at -20°C till use.

- **Preparation of the conjugate and substrate:** do not dilute. Ready to use solution, add 100 µl/well

VI. PREPARATION OF SAMPLES:

For diagnosis purposes, samples must be whole blood, spleen or cellular culture.

For whole blood (with heparin or EDTA), before use, samples must be frozen and thawed once to avoid false positives results.

Samples are used at 1/2 dilution (50 µl of diluent and 50µl of sample/well).

For spleen, a 1/10 dilution must be done (w/v) using the diluent supplied (i.e. 1g of sample in 10 ml of diluent). Make a good homogenisation of the sample dilution using a stomacher or any other mechanic homogenizer. Filter or centrifuge the homogenized material to remove tissue debris and add 100 µl to the well

VII. TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. Add 100 µl of positive control to two wells (e.g. A1 and B1), and 100 µl of the negative control (e.g. A2 and B2). Add 100 µl of samples to test on each remainder wells (if the sample is whole blood you can make the dilution directly in the well adding in first place the diluent and in second place the sample). We recommend the use of two wells per sample. Seal the plate and incubate **for 1 h at 37°C**.
3. **Empty the wells into a receptacle containing 0.1 M NaOH** and wash four times as specified in point IV.
4. Add 100 µl of specific conjugate to each well. Seal the plate and incubate for **1 hour at room temperature (20-25°C)**.
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate to each well. Keep the plate for **15 min at room temperature**.
7. Add 100 µl of stop solution to each well.
8. Read the OD of each well at 405nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

Validation of the test

The test could be considered valid when:

The OD of the Positive control is higher that 1.0

The OD of the Negative control is lower that 0,150

Interpretation:

When you are running duplicates, OD of the sample, will be calculated as the arithmetic mean of OD values in both wells.

- Samples with an OD equal or higher than 0.250 must be considered as **POSITIVES** to ASFV.
- Samples with OD equal or lower that 0,200 must be considered as **NEGATIVE** to ASFV.
- Samples with OD values between both values are considered as **DOUBTFUL**. We recommend re-testing these animals one more time or applying a different technique to check this sample

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA
APLICADA, S.A.
Av. Institución Libre de
Enseñanza, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by: