
INgezim® FPV IC

R.16.FPV.K42

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

Ed.110418

INgezim® FPV IC R.16.FPV.K42

INgezim® FPV IC is a direct immunochromatographic assay, based on Lateral Flow technique, which uses 2 monoclonal antibodies (MAbs) specific of Feline Panleukopenia virus (FPV).

TECHNICAL BASIS OF THE KIT

The device consists on a plastic casing with 2 windows:

- **Window for sample addition:** Contains a MAb specific FPV and a Control Protein conjugated with colloidal gold.
- **Window for results interpretation:** Contains a test line (T) with a FPV specific MAb and a control line (C) with a MAb specific to the control protein.

APPLICATION

Detection of Feline Panleukopenia virus in biological samples (faeces and intestinal fluid) of cats.

INTERPRETATION OF THE RESULTS

Positive sample: Pink line in the Test zone and a pink line in the Control zone.

Negative sample: Only a pink line in the Control zone.

Not valid assay: No pink line in the Control zone

SENSITIVITY AND SPECIFICITY

Correspondence with PCR technique.

154 faeces samples of cats collected from the clinicians in Italy. The results were compared with the ones obtained by PCR showing 98.2% and 100% sensitivity and specificity respectively regarding PCR.

		PCR		
		POS	NEG	TOT
FPV IC	POS	55	0	55
	NEG	1	98	99
	TOT	56	98	154

COMPOSITION OF THE KIT

- Diagnostic devices individually wrapped
- Dropper vial with dilution buffer.

PRODUCT MANUFACTURED BY
AGROLABO AND DISTRIBUTED BY
INGENASA

SHELF LIFE: 18 months
Stored: 4°C-25°C