INgezim Brucella IC

R.15.BC.K41

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
INgezim Brucella IC

R.15.BC.K41

INgezim Brucella IC is an indirect immunochromatographic assay, based on Lateral Flow technique, which uses a monoclonal antibody specific of canine IgG and a specific antigen (Ag) of *Brucella canis*.

TECHNICAL BASIS OF THE KIT

The device consists on a plastic casing with 2 windows:

- **Window for sample addition:** Contains a MAb specific of canine IgG and a Control Protein. All of them labelled with colloidal gold.
- **Window for results interpretation:** Contains a test line (T) with a specific antigen of *B. canis* and a control line (C) with a MAb specific to the control protein

APPLICATION

Detection of antibodies specific of *Brucella canis* in serum, plasma and blood samples of dogs.

INTERPRETATION OF THE RESULTS

Positive sample: Pink line in the Test zone and a pink line in the Control zone.

Negative sample: Only a pink line in the Control zone.

Not valid assay: No pink line in the Control zone

SENSITIVITY AND SPECIFICITY

Correspondence with culture.

176 canine serum and plasma samples from different clinicians were analyzed. Results were compared with the ones obtained by culture isolation. The assay shows 94% and 100% sensitivity and specificity respectively.

		CULTURE		
		POS	NEG	TOT
EHRLI CROM	POS	48	0	48
	NEG	3	125	128
	TOT	51	125	176

COMPOSITION OF THE KIT

- Diagnostic devices individually wrapped
- Diluent tubes to prepare the dilution of the faeces.
- Plastic pipettes

PRODUCT MANUFACTURED BY
AGROLABO AND DISTRIBUTED BY
INGENASA

SHELF LIFE: 18 months
Stored: 4°C-25°C

Ed.110418