INgezim® CSF Compac

R.11.CSF.K3

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INgezim® CSF Compac

R.11.CSF.K3

INgezim® CSF Compac is an enzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific to Classical Swine Fever Virus (CSFV) E2 protein as conjugate and a recombinant E2 protein as antigen coating the plate.


KIT FEATURES

APPLICATION

Detection of specific antibodies to Classical Swine Fever Virus, in porcine sera samples.

TECHNICAL BASE

- Plates are supplied coated with CSFV antigen (E2 recombinant protein). Samples are added to the wells and incubated.
- If the samples contain specific antibodies to CSFV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific of E2 protein, conjugated with peroxidase, AcM-PO) is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected or vaccinated animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

RESULTS INTERPRETATION

The assay establishes one cut off, which will classify the samples as *Positive* or *Negative*, depending on the value of the blocking percentage of the sample in the assay.

ASSAY VALIDATION

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

932 sera were analysed: 902 from Spanish CSF free areas and a reference panel classified by Seroneutralization (SN) (Chapter VII of Commission Decision 2002/106/EC): 11 early phase infected (≤ 21 d.p.i), 13 convalescent phase (> 21 d.p.i) and 6 sera positive to Pestiviruses different than CSFV. The assay showed 95% sensitivity and 99.9% specificity.

ANALYTICAL SENSITIVITY AND SPECIFICITY

30 sera belonging to the reference panel described below were analysed (Chapter VII of the Commission Decision 2002/106/EC): 11 early phase infected (≤ 21 d.p.i), 13 convalescent phase (> 21 d.p.i) and 6 sera positive to Pestiviruses different than CSFV. Results obtained indicate that the assay is able to detect sera with SN titers of 1/30 and does not show cross reaction with antibodies specific to other pestiviruses like BDV and BVDV.

COMPARATIVE WITH OTHER ASSAYS AVAILABLE

• Study 1:

INgezim® CSF COMPAC has been compared with 2 assays available in the market (PrioCheck CSFV 2.0 Ab y HerdCheck CSFV Ab) regarding to their behaviour with respect to SN. For that propose, the panel described above, was used. Results obtained are shown below:

CORRESPONDENCE WITH SN			
	PrioCheck CSFV Ab 2.0	HerdCheck CSFV Ab	INgezim® CSF Compac
Specificity	100%	100%	100%
Sensitivity	95.4%	95.8%	95.8%

• Study 2:

INgezim® CSF COMPAC was compared with 2 assays available in the market using a panel of vaccinated animals. Results indicated that INgezim® CSF COMPAC showed 100% sensitivity and specificity. In this study it was also confirmed that the assay detects specific antibodies at 4 weeks after vaccination in 62.5% of cases. The table below indicates the comparative with the other 2 assays.

	INGENASA		TEST I		TEST II	
	POS	NEG	POS	NEG	POS	NEG
Vaccinated and revaccinated	64	0	62	2	8	56
Vaccinated	0	40	2	38	2	38
4 weeks post vaccination	50	30	50	30	22	58
Total	114	70	114	70	32	152

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

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Eurofins-INGENASA, S.A.