INgezim® PPA Compac

R.11.PPA.K3

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INgezim® PPA Compac

R.11.PPA.K3

INgezim® PPA Compac is an enzymatic assay based on a blocking ELISA assay technique, which uses a monoclonal antibody (MAb) specific to African Swine Fever Virus (ASFV) VP72 protein.


KIT FEATURES

APPLICATION

Detection of specific antibodies to African Swine Fever Virus, in porcine sera samples.

TECHNICAL BASE

- Plates are supplied coated with ASFV antigen (VP72 protein). Samples are added to the wells and incubated.
- If the samples contain specific antibodies to ASFV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific of VP72 protein, conjugated with peroxidase, AcM-PO) is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the protein. In case the sample contains antibodies blocking the antigen (infected or vaccinated animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

RESULTS INTERPRETATION

The assay establishes two cut offs, which will classify the samples as **Positive** or **Negative**, depending on the value of the optical density of the sample in the assay.

ASSAY VALIDATION

DIAGNOSTIC SENSITIVITY


208 sera from endemic countries in West and East Africa and from old outbreaks in Spain were analysed. These sera had previously been classified by the OIE Reference Indirect ELISA as positive samples. The results obtained indicated a 99% sensitivity respect the OIE Technique.

DIAGNOSTIC SPECIFICITY

- Testing sera from endemic areas in Africa.** 167 African pig sera from ASFV endemic areas were analysed. These sera had been classified as negative by the OIE reference Indirect ELISA. The results obtained showed 100% specificity.
- Testing sera from ASFV non infected areas.** 1043 sera from ASFV non-infected areas were analysed. These sera had been classified as negative by the OIE reference Indirect ELISA. The results obtained showed 100% specificity.
- Testing problematic sera.** A set of 23 samples with negative results by Immunoblotting and doubtful results by the Indirect ELISA was analysed. The results obtained by INgezim® PPA COMPAC showed a correspondence of 100% with the Immunoblotting technique.

ANALYTICAL SENSITIVITY

- Study 1.** Samples from pigs experimentally infected with different isolates of ASFV were analyzed. The isolates used for the experimental infection were: E70/E75 p72 genotype I (Spain); p72 genotype I (Sardinia); Ken05.Tk1 p72 genotype X (Kenya); p72 genotype IX (Kenya); p72 genotype II (Armenia); Malau 83, genotype VIII and Rep. Dominicana 78, genotype I. **The results obtained indicated that the assay is able to detect antibodies specific of all of these isolates. These antibodies are detected between days 10 and 25 post infection, depending on the isolate used.**
- Study 2.** Different extractions (d.p.i. 0 to 59) from 30 experimentally infected animals (Benin strain, genotype I) were analyzed using INgezim® PPA Compac. Results were compared with 2 indirect ELISAs for detection of ASFV specific IgM and IgG respectively. This study allows to observe the good performance of the assay for detection of IgM, confirming its precocity in detection of ASFV specific antibodies (since day 10 p.i.).

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

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