INgezim® PARVO CANINO

R.15.CPV.K1

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INgezim® PARVO CANINO

R.15.CPV.K1

INgezim® PARVO CANINO is based on an indirect enzymatic immunoassay (indirect ELISA), which uses a monoclonal antibody (MAb), specific to canine immunoglobulins (Igs), and a canine Parvovirus (CPV) VP2 recombinant protein as antigen.


KIT FEATURES

APPLICATION

Detection and/or titration of specific antibodies to CPV in canine sera samples.

TECHNICAL BASE

- Plates are supplied coated with CPV antigen. Serum samples are added and incubated.
- If the samples contain specific antibodies to CPV, they will bind to the antigen.
- When the conjugate (MAb-PO specific to canine immunoglobulins) is added, it will bind to the Igs previously bound to the antigen. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.


RESULTS INTERPRETATION

The assay establishes a cut off, which will classify the samples as *Positive* or *Negative*, depending on the Positive Index (PI) of the sample in the assay. The titre of each sample will be the last dilution at which the IP value is higher than the cut off.

ASSAY VALIDATION

Analytical sensitivity in relation to IHA:

- Study 1:** 20 sera of animals experimentally vaccinated and challenged were analysed. 19 extractions were taken and analysed by both assays. The results obtained indicated that ELISA is able to detect specific antibodies of CPV at day 10p.v. in 35% of cases and IHA in 25%.
- Study 2:** 10 sera of experimentally vaccinated dogs were analysed. 20 extractions were taken and analysed by both assays. The results obtained indicated that the ELISA is able to detect specific antibodies of CPV at day 10p.v. in 20% of cases and IHA 0%.


Correlation with inhibition haemagglutination test (IHA). Field sera


100 canine field sera were analysed in this study. The results obtained indicated 95% sensitivity and a specificity percentage higher than 99%.

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- Field studies:** The graph indicates the correlation between titres of both techniques in 61 animals analysed.

- Experimental studies:** The graph indicates the titre correlation in a study made with 20 dogs experimentally vaccinated and challenged at day 73 post vaccination.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

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