INgezim® BRUCELLA Compac 2.0

R.10.BRU.K3

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Referencias: García-Bocanegra et al. Vet. Evaluation of different enzyme-linked immunosorbent assays for the diagnosis of brucellosis due to *Brucella melitensis* in sheep J. 2014 Mar; 199(3):439-45.

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INgezim® BRUCELLA Compac 2.0

R.10.BRU.K3

INgezim® BRUCELLA Compac is a multispecies enzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (Mab) specific to the epitope C of LPS of Brucella.


KIT FEATURES

APPLICATION

Detection of specific antibodies to LPS of *Brucella spp.*, in bovine, ovine, caprine and porcine sera samples.

TECHNICAL BASE

- Plates are supplied coated with LPS antigen of *Brucella abortus*. Serum samples are added to the wells and incubated.
- If the samples contain specific antibodies to *Brucella spp.*, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific to LPS, conjugated with peroxidase, AcM-PO) is added, it will bind to the protein only if there are no antibodies in the sample blocking the antigen (negative animals). In case there were antibodies blocking the antigen (infected or vaccinated animals), the conjugate would not be able to bind it. This binding is revealed by colorimetric reaction after addition of the substrate.

RESULTS INTERPRETATION

The assay establishes a cut off, which will classify the samples as **Positive** or **Negative**, depending on the value of inhibition percentage (IP) of the sample in the assay.

ASSAY VALIDATION

ANALYTICAL SENSITIVITY

The assay fulfils the OIE requirements and the European Directives regarding bovine, small ruminant and porcine International Reference Sera.

DIAGNOSTIC SENSITIVITY

- 252 sera (74 bovine, 66 ovine, 46 caprine and 66 porcine), positive by cell culture, were analysed.
- 184 sera, positive by other techniques (124 bovine, 24 ovine, 22 goats, 14 porcine).

The diagnostic sensitivity was 98% for cattle and 99% for sheep, goats and swine.

RING TEST

This assay has been evaluated with periodical Ring Tests organized by the Animal Health Laboratory in Maisons-Alfort, France (ANSES) in years 2012, 2014, 2016 and 2018. The results indicated 100% correspondence with the expected results.

References: García-Bocanegra et al. Vet. Evaluation of different enzyme-linked immunosorbent assays for the diagnosis of brucellosis due to *Brucella melitensis* in sheep J. 2014 Mar; 199(3):439-45.

ANALYTICAL SPECIFICITY

100% using negative reference sera, 97% analysing animals infected with *Yersinia enterocolitica* and 95.4% for RB51 positive culture aborted animals.

DIAGNOSTIC SPECIFICITY

- 5278 sera (2817 bovine, 810 ovine, 719 goats and 932 porcine) from Brucella free herds were analysed. The diagnostic specificity was 99.9%.
- 197 sera (bovine, goats, ovine and porcine) culture negative were analysed. The diagnostic specificity was 99.9%.

STUDIES IN CAMELIDS

63 camelid sera (5 camel & 58 alpacas) have been studied. Results showed 100% correspondence with the expected.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

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