INgezim® COVID 19 S VET

R.10.CoS.K1

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INgezim® COVID 19 S VET

R.10.CoS.K1

INgezim® COVID 19 S VET is based on an indirect ELISA assay that uses SARS-CoV-2 S protein as antigen and recombinant protein A/G conjugated with peroxidase (PO) as conjugate.


KIT FEATURES

APPLICATION

Detection of specific antibodies against S protein of the SARS-CoV-2 virus in serum and plasma samples from minks, ferrets, cats and dogs.

TECHNICAL BASE

1. Plates are supplied coated with antigen (SARS-CoV-2 recombinant S protein). Samples are added to the wells and incubated.
2. If the samples contain specific antibodies to SARS-CoV-2 S protein, they will bind to the antigen.
3. When the conjugate (recombinant protein A/G, conjugated with peroxidase) is added, it will bind to the antibodies bound to the antigen. This binding is revealed by colorimetric reaction after addition of the substrate.

RESULTS INTERPRETATION

A cut off has been established for each animal species, which will classify the samples as **Positive** or **Negative**, depending on the value of optical density of the sample in the assay.

ASSAY VALIDATION

ANALYTICAL SPECIFICITY

- 70 **dog sera** positive for antibodies to canine coronavirus were analysed. The obtained results indicated that **there is no cross-reaction**.
- 9 **cat sera** positive for antibodies to feline coronavirus were analysed. The obtained results indicated that **there is no cross-reaction**.

ANALYTICAL SENSITIVITY

Two sera from experimentally infected cats were analysed, which were sampled at different days post infection. The obtained results indicated that the assay can detect specific antibodies from day 10 post infection.

DIAGNOSTIC SPECIFICITY & SENSITIVITY

In order to determine the diagnostic specificity and sensitivity in different types of samples, the following sera were analysed:

- Mustelids samples:
 - 243 negative mink sera and 31 negative ferret sera
 - 211 positive mink sera
- Dog samples:
 - 355 negative sera
 - 7 positive sera
- Cat samples:
 - 223 negative sera
 - 11 positive sera

	Sensitivity	Specificity
Mustelids samples	98,1%	>99,9%
Dog samples	>99,9%	98,9%
Cat samples	>99,9%	99,6%

CORRESPONDENCE WITH SERONEUTRALIZATION

- 21 **dog samples** previously characterized by seroneutralization were analysed. The results showed **95% correspondence** between both techniques.
- 17 **cat samples** previously characterized by seroneutralization were analysed. The results showed higher than **99.9% correspondence** between both techniques.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

EXPIRATION: 6 MONTHS (*). Stored at 2°C-8°C
 (*) Expiration studies in progress.

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