INgezim® CORONA DIFERENCIAL 2.0

R.11.DIF.K3

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INgezim® CORONA DIFERENCIAL 2.0

R.11.DIF.K3

INgezim® CORONA DIFERENCIAL 2.0 is based on an enzymatic immunoassay (Capture Blocking ELISA) which uses as conjugates two specific monoclonal antibodies (MAB); one specific to site B of coronavirus S protein (generic epitope of coronavirus) and the other one specific to site Ac of TGEV S protein (specific epitope of TGEV). A third MAB is used as a capture reagent of the S protein.


KIT FEATURES

APPLICATION

Detection and differentiation of Transmissible Gastroenteritis Virus (TGEV) and Porcine Respiratory Coronavirus (PRCV) specific antibodies in porcine sera samples.

TECHNICAL BASE

- Plates are coated with recombinant S protein of TGEV which is captured by a specific monoclonal antibody. Each serum sample must be added to two wells.
- If the serum sample contains specific antibodies to the Porcine Respiratory Coronavirus (PRCV), they will only bind to the generic epitopes, but if it contains specific antibodies to the Transmissible Gastroenteritis Virus (TGEV), they will bind to all of the epitopes present.
- After an incubation period, a coronavirus specific monoclonal antibody-HRPO (Conjugate A) must be added to one well and a specific TGEV MAb-HRPO (Conjugate B) to the other one. They will only bind to free epitopes. The presence or absence of the labelled MAB can be detected by adding a substrate which, in presence of the peroxidase, will develop a colorimetric reaction.

RESULTS INTERPRETATION

Three cut off are used for the results interpretation: positive/negative to coronavirus, positive to TGEV and positive to PRCV/negative to TGEV. The combination of these three Cut Offs allows differentiating between TGEV positive samples, PRCV positive samples and Coronavirus negative samples.

ASSAY VALIDATION

EXTERNAL EVALUATION

A study with 411 samples, divided into 4 groups, was made:

- 94 samples from 16 herds serologically negative to PRCV and TGEV.
- 124 samples from 21 herds serologically positive to PRCV.
- 122 samples from 17 herds serologically positive to TGEV
- 71 samples from 22 diagnostic cases

The results obtained indicated **94% sensitivity** and **98.2% specificity**.

INTERNAL EVALUATION

Comparison with an ELISA that uses complete TGEV virus as antigen.

A set of 665 sera were analysed and the results obtained in both assays were compared:

		INGEZIM CORONA DIFERENCIAL 2.0				
		PRVC	TGEV	NEGATIVE	DOUBTFUL	TOTAL
ELISA	PRCV	498			5	503
	TGEV	1	31		1	33
	NEGATIVE	11		112		123
	DOUBTFUL	6				6
	TOTAL	516	31	112	6	665

Those samples with discrepant results showed OD values close to the cut off in the ELISA that used complete virus as antigen. The use of recombinant protein as antigen allows the detection of specific antibodies to TGEV, avoiding false positive results.

COMPOSICION DEL KIT

- Microtitration plates of 96 wells.
- Vials with Positive Control for PRCV
- Vials with Positive Control for TGEV
- Vials with Negative Control for Porcine Coronavirus
- Vials with Peroxidase Conjugate A (TGEV+PRCV)
- Vials with Peroxidase Conjugate B (TGEV)
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

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