INgezim® BLV Confirmation

R.12.BLV.K1

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INgezim® BLV Confirmation

R.12.BLV.K1

INgezim® BLV Confirmation is an immunoenzymatic assay based on an indirect ELISA technique, which uses a monoclonal antibody (MAb) specific of bovine immunoglobulins. This assay uses plates coated with positive antigen (inactivated BLV) and negative antigen.


KIT FEATURES

APPLICATION

Confirmatory assay for detection of specific antibodies to Bovine Leukaemia Virus (BLV) in bovine serum samples.

TECHNICAL BASE

- Plates are supplied coated with BLV antigen and negative antigen, respectively. Samples are added to the wells and incubated.
- If the samples contain specific antibodies to BLV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific of bovine immunoglobulins, conjugated with peroxidase, AcM-PO) is added, it will bind to the IgGs bound to the antigen. This binding is detected by the development of a colorimetric reaction after the addition of the substrate.

RESULTS INTERPRETATION

The assay establishes one cut off, which will classify the samples as *Positive* or *Negative*, depending on the value of the optical density of the sample in the assay.

ASSAY VALIDATION

USE OF FIELD SERA

A study was made using a set of 400 sera previously classified by INgezim® BLV COMPAC:

- 29 sera positive from Canada and Europe
- 371 sera negative from Holland, Spain, Poland, and Canada

Classified	Nº	POS	DUD	NEG
Negatives	371	4	0	367
Positives	29	28	0	1

Correspondence between both assays was 98.7%.

USE OF O.I.E. REFERENCE SERA

In order to verify the sensitivity of the assay, E-4 and E05 Community Reference Sera were requested, and, following the O.I.E. legislation, they were assayed at a 1/10 dilution in negative serum (as weak positive serum) and diluted 1/100 in negative sera. Positive results were obtained in all cases. The OIE Negative Reference Serum showed negative results in the same assay. Therefore, we conclude that INGEZIM BLV Confirmation maintains the level of sensitivity and specificity required by the O.I.E legislation, and therefore can be used as a confirmation assay.

KIT COMPOSITION

- Microtitration plates coated with BLV antigen
- Microtitration plates coated with negative antigen
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

Spanish registration nº 3471RD

EXPIRATION: 18 MONTHS. Stored at 2°C-8°C

Eurofins-INGENASA, S.A
 Avda. de la Institución Libre de Enseñanza 39, 8º
 28037 MADRID (SPAIN)
 Phone: (+34)91 3680501
www.ingenasa.com