# INgezim® PPA DAS 2.0

R.11.PPA.K2

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# INgezim® PPA DAS 2.0

R.11.PPA.K2

INgezim® PPA DAS 2.0 is an enzymatic assay based on the Double Antibody ELISA technique, which uses monoclonal antibodies (Mab) specific to the African Swine Fever Virus (ASFV) VP72 protein.


## KIT FEATURES

### APPLICATION

Detection of African Swine Fever Virus (ASFV) in whole blood, cellular culture and spleen samples.

### TECHNICAL BASE

1. Plates are supplied coated with antibody (monoclonal antibody specific to the VP72 protein of ASFV). Samples are added to the wells and incubated.
2. If the samples contain antigens, they will bind to the specific antibodies for the VP72 protein which are coated in the plates.
3. When the conjugate (monoclonal antibody specific to a different epitope of VP72, conjugated with peroxidase, AcM-PO) is added, it will bind to the antigen bound to the antibodies coated to the plate. This binding is revealed by colorimetric reaction after addition of the substrate.

### RESULTS INTERPRETATION

The assay establishes two cut offs, which will classify the samples as *Positive* or *Negative*, depending on the value of the optical density of the sample in the assay, considering a range of optical densities close to the cut offs as *Doubtful* results.

## ASSAY VALIDATION

### DETECTION LIMIT


Different amounts of u.f.p./ml (plate forming units/ml) OF ASFV inoculum were titrated. The obtained results indicated that the assay is capable of detecting 5x10<sup>4</sup> u.f.p./ml of ASFV.

### SPLEEN SAMPLES

- **Diagnostic Specificity.** A panel of spleen samples from uninfected animals was analysed. The assay classified 100% of these samples as negative.
- **Diagnostic Sensitivity.** A panel of spleen samples from infected animals was analysed. The assay classified 100% of these samples as positive.

### WHOLE BLOOD SAMPLES

- **Analytical Sensitivity.** Whole blood samples from 4 experimentally inoculated animals were analysed. The obtained results indicated that the assay can detect ASFV antigen from day 7 post-infection.
- **Detection of different isolates.** Samples of animals infected with different ASFV isolates belonging to genotypes I, II, IX & X were studied. The assay can detect all isolates analysed.
- **Diagnostic Specificity.**
  - **Internal study.** 278 negative samples from Spain (ASFV free country) were analysed. All of them were also analysed by INgezim® ASFV CROM Ag, confirming their negative status. The obtained results indicated 99.6% specificity.
  - **External study.** 33 blood samples negative by PCR-UPL and isolation were analysed. The obtained results indicated 97% specificity with respect to PCR-UPL and isolation in this study.
- **Diagnostic Sensitivity. External study.** The EU Reference Laboratory (CISA-INIA) carried out a study using 146 blood samples from animals experimentally infected with different isolates and positive by UPL-PCR. The sensitivity with respect to UPL-PCR depends on the degree of positivity of the sample. For Ct < 20%, 98% sensitivity was obtained with respect to UPL-PCR, for Ct ≥ 20 < 25, 62% sensitivity was obtained with respect to UPL-PCR.

## KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (ABTS)
- Bottles with Stop Solution

Spanish registration n° 3465RD

EXPIRATION: 18 MONTHS. Stored at 2°C-8°C

**Eurofins-INGENASA**  
 Avda. de la Institución Libre de Enseñanza 39, 8º  
 28037 MADRID (SPAIN)  
 Phone: (+34)91 3680501  
[www.ingenasa.com](http://www.ingenasa.com)

IT-73840  
IT-73780

9191.INGE

9175.ING2