INgezim® BLV Compac 2.0

R.12.BLV.K3

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INgezim® BLV Compac 2.0

R.12.BLV.K3

INgezim® BLV Compac 2.0 is an immunoenzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific for Bovine Leukaemia (BLV) gp51 protein.


KIT FEATURES

APPLICATION

Detection of specific antibodies to BLV gp51 protein in serum samples from cattle individually assayed or pooled up to 10 and milk samples individually assayed.

TECHNICAL BASE

- Plates are supplied coated with BLV antigen (inactivated BLV gp51 protein). Samples are added to the wells and incubated.
- If the samples contain specific antibodies to gp51 protein of BLV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific of gp51 protein, conjugated with peroxidase, AcM-PO) is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the antigen. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

RESULTS INTERPRETATION

The assay establishes two cut offs, which will classify the samples as *Positive* or *Negative*, depending on the value of the optical density of the sample in the assay, considering a range of optical densities close to the cut offs as *Doubtful* results.

ASSAY VALIDATION

SENSITIVITY

Using O.I.E. reference sera

O.I.E. reference sera E-4 and E05 were assayed without prior dilution (as strong positive serum), at 1/10 dilution in negative serum (as weak positive serum) and at dilution 1/250 in negative milk (as individual milk sample). Positive results were obtained in all cases, concluding that INgezim® BLV Compac 2.0 maintains the level of sensitivity required by the O.I.E. (Manual Standards for Diagnosis Test and Vaccines 6th Edition 2008) for individual and pools of 10 sera and for individual milk.

SPECIFICITY

In order to determine the specificity of the assay, a set of 354 samples of bovine sera from four different herds, were analysed. The results obtained indicated **>99% specificity**.

FIELD TRIALS

Field sera were analysed comparing the results obtained between INgezim® BLV Compac 2.0 and the reference technique, AGID. The results showed **>99% correlation**.

EXTERNAL VALIDATION

The assay has been validated by the OIE Reference Laboratory for BLV in Leipzig (Dr. Thomas W. Vahlenkamp) (2016)

The assay was validated by the National Veterinary Research Institute (PIWet, Poland) in 2005.

The assay has been registered in Germany with n° FLI-C 033.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

Spanish registration n° 808RD

EXPIRATION: 18 MONTHS. Stored at 2°C-8°C

Eurofins-INGENASA, S.A
 Avda. de la Institución Libre de Enseñanza 39, 8º
 28037 MADRID (SPAIN)
 Phone: (+34)91 3680501
www.ingenasa.com

IT-73840
IT-73780

9191.INGE 9175.ING2