## INgezim® *Brucella* Small Ruminants

R.13.BM.K1

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## INgezim® *Brucella* Small Ruminants

R.13.BM.K1

INgezim® *Brucella* Small Ruminants is an immunoenzymatic assay based on an indirect ELISA technique, which uses a monoclonal antibody (MAb) specific for ruminant immunoglobulins and *Brucella mellitensis* (LPS) antigen.


### KIT FEATURES

#### APPLICATION

Detection of specific antibodies to *Brucella* in small ruminants serum samples individually assayed.

#### TECHNICAL BASE

1. Plates are supplied coated with inactivated *Brucella mellitensis* antigen (LPS). Samples are added to the wells and incubated.
2. If the samples contain specific antibodies to *Brucella*, they will bind to the antigen.
3. When the conjugate (monoclonal antibody specific to ruminant immunoglobulins, conjugated with peroxidase, AcM-PO) is added, it will bind to the IgGs bound to the antigen. This binding is detected by the development of a colorimetric reaction after the addition of the substrate.

#### RESULTS INTERPRETATION

The assay establishes one cut off, which will classify the samples as **Positive** or **Negative**, depending on the value of the optical density of the sample in the assay.

### ASSAY VALIDATION

#### Correspondence with Rose of Bengal (RB) and Complement Fixation (CF)

In order to determine the correspondence with both reference techniques, a set of 345 sera from animals (goats and sheep) under eradication programmes were analysed. The results obtained indicated 96.2% correspondence with CF and 95% with RB.

#### Correspondence with culture

- **Study 1.** 127 samples from non-vaccinated sheep in infected herds were analysed. These samples were positive by culture. The results obtained indicated that the sensitivity of the assay was 96.1% being this percentage higher than the one obtained by RB and CF.
- **Study 2.** 275 samples from sheep and goats, negative and positive by isolation were analysed. The results obtained indicated 98.9% sensitivity for sheep and 97.7% for goats, and >99% specificity.

#### Analytical sensitivity and specificity. Use of OIE Reference Sera

This assay has been validated with the new International OIE Reference Sera EUSBSS and EUGBSS, which derive from the ISaBmS serum (against which the assay has also been validated).

EUGBSS and EUSBSS must be positive at least at dilution 1/8 and negative at dilution 1/256 as required by the OIE. The assay detects as positive both sera at dilution 1/64, being dilution 1/250 negative. Regarding the International Standard ISaBmS, the assay detects as positive the dilution 1/250 and as negative the dilution 1/500.

#### Specificity of the assay

480 sera from 6 *Brucella* free herds were analysed. The results obtained indicated >99.9% specificity.

### KIT COMPOSITION

- 96 wells microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

Spanish registration nº 3057RD

EXPIRATION: 24 MONTHS. Stored at 2°C-8°C

**Eurofins-INGENASA, S.A**  
 Avda. de la Institución Libre de Enseñanza 39, 8º  
 28037 MADRID (SPAIN)  
 Phone: (+34)91 3680501  
[www.ingenasa.com](http://www.ingenasa.com)

IT-73840  
IT-73780

9191.INGE

9175.ING2