INgezim® PPR Compac

R.13.PPR.K3

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INgezim® PPR Compac

R.13.PPR.K3

INgezim® PPR Compac is an immunoenzymatic assay based on a blocking ELISA technique, which uses a monoclonal antibody (MAb) specific to Pest des Petites Ruminant Virus (PPRV) N protein and the recombinant N protein of the PPRV.


KIT FEATURES

APPLICATION

Detection of specific antibodies to the N protein of PPRV. The detection of antibodies does not depend on the animal species and can therefore be used with sera samples of all species affected by PPRV.

TECHNICAL BASE

- Plates are supplied coated with antigen (recombinant N protein of PPRV). Samples are added to the wells and incubated.
- If the samples contain specific antibodies to PPRV, they will bind to the antigen.
- When the conjugate (monoclonal antibody specific to PPRV N protein, conjugated with peroxidase, AcM-PO) is added, only if there are no antibodies in the sample blocking the antigen (negative animals), it will bind to the antigen. In case the sample contains antibodies blocking the antigen (infected animals), the conjugate will not be able to bind to it. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

RESULTS INTERPRETATION

The assay establishes one cut off, which will classify the samples as *Positive* or *Negative*, depending on the Inhibition Percentage (IP) value of the sample in the assay.

ASSAY VALIDATION

DIAGNOSTIC SPECIFICITY

2698 animals (sheep & goats) from different PPRV free areas were analysed. Results obtained indicated **97.5% specificity**.

DIAGNOSTIC SENSITIVITY

113 sera from different sources, both field and experimentally vaccinated ones, were analysed. Field sera were previously classified by ID Screen® PPR Competition. Results obtained indicated **88% relative sensitivity**.



	PPR positive sera	Positives by INgezim PPR COMPAC
Cameron (field)	46	40
Morocco (field)	53	47
Spain (experimental)	14	12
Total	113	99

ANALYTICAL SENSITIVITY

Sera of experimentally vaccinated and infected animals were analysed:

- 7 one year old sheep, subcutaneously vaccinated with attenuated vaccine Nigeria75/1 (linage II). Extractions were made at days 4 and 7 p.v.
- 1 sheep one year old, intravenously inoculated with the field virulent strain Morocco 2008 (linage IV). Extractions at days 0, 2, 4, 7, 10, 14, 17, 21, 24, 28 & 35 p.i.23 were made.

The samples were analysed by INgezim® PPR Compac. Results obtained indicated that the assay detects antibodies specific to PPRV at days 7 p.v. and p.i. in 85% of cases.

KIT COMPOSITION

- 96 well microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with Conjugate
- Bottles with Washing Solution
- Bottles with Diluent
- Bottles with substrate (TMB)
- Bottles with Stop Solution

Spanish registration nº 3780RD

EXPIRATION: 24 MONTHS. Stored at 2°C-8°C

Eurofins-INGENASA, S.A
 Avda. de la Institución Libre de Enseñanza 39, 8º
 28037 MADRID (SPAIN)
 Phone: (+34)91 3680501
www.ingenasa.com

IT-73840
IT-73780

9191.INGE 9175.ING2