

INgezim BESNOITIA

R.12.BES.K1

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

INgezim BESNOITIA

R.12.BES.K1

INgezim BESNOITIA is an indirect immunoenzymatic assay that uses an extract of proteins from the tachyzoite stage of *Besnoitia* to coat the plate and a monoclonal antibody specific to bovine IgG as conjugate.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with an extract of proteins of *Besnoitia*. Serum samples are added to the wells and incubated.
2. If the sample has specific antibodies to *Besnoitia*, they will bind to the antigen.
3. When a HRPO conjugated is added, it will bind to the antibodies present in the sample (positive serum). This binding is detected by adding a specific substrate.

APPLICATION

Detection of antibodies specific of *Besnoitia* in Bovine serum.

INTERPRETATION OF THE RESULTS

The assay establishes two cut offs. Samples will be considered **Positive** if its IRPC value is higher than the positive cut off; **Negative**, if its IRPC value is lower than the negative cut off. IRPC values between both cut offs will be considered **Doubtful**.

SENSITIVITY

88 sera from a herd suspected of being infected with *Besnoitia* were analyzed. The assay was able to detect 80% of these animals as positive.

CORRESPONDENCE WITH OTHER ASSAYS

- **STUDY 3:** the results obtained for 241 sera which were analyzed by INGEZIM *Besnoitia*, Indirect ELISA^a, Westernblot (WB) to Tachyzoite and to a Bradizoite, were compared. The correspondence between INGEZIM and the other assays is shown below:

	WB TACHYZOITE	WB BRADIZOITE	INDIRECT ELISA ^a
Sensitivity	97%	95%	96%
Specificity	83%	89%	92%
Correspondence	89%	91%	94%

^a Fernández-García et al. Veterinary Record (2010); 166: 818-822

- **STUDY 2:** 295 field samples were analyzed by two different assays: INGEZIM *Besnoitia* and PrioCHECK *Besnoitia*. These samples were also checked for *Toxoplasma* and *Neospora*, being 251 positive or border to *Toxoplasma* and 4 to *Neospora*. The results obtained are shown below:

	INGEZIM BESNOITIA				
	Neg.	Pos.	Doubt.	Total	
PRIOCHECK <i>Besnoitia</i>	Neg.	103	4	0	107
	Pos.	20	65	3	88
	Doubt.	83	12	5	100
	Total	206	81	8	295

*72% of toxoplasma positive were positive to *Besnoitia* by PrioCheck but only 31% were positive by INGEZIM.

SPECIFICITY

145 sera of animals negative to *Besnoitia* and positive to *Neospora* (natural or experimental infection) were analyzed. Results indicated 98% specificity.

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate buffer
- Bottle with substrate (ABTS)

PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **12 months**
Stored at 2°C-8°C

Ed.020217