# VetLine Bovine Haptoglobin

**ELISA**

**RUO**

Enzyme immunoassay for the quantitative determination of Haptoglobin in bovine serum or milk.

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---

Product Number:      BHAPV4040 (96 Determinations)

---

## ENGLISH

### 1. INTRODUCTION

---

Haptoglobin (Hp) is the most studied acute phase protein in cattle. Numerous studies have shown that the concentration of protein in the serum during infectious and inflammatory diseases increases sharply. During mastitis, the Hp concentration in milk is also greatly increased. In the course of peripartum immunosuppression and early lactation, dairy cows are also particularly susceptible to infectious diseases such as uterine, claw or respiratory infections. Not infrequently, a relocation of the abomasum occurs. It could be shown that even during these diseases the haptoglobin concentration in the milk is increased.

The determination of Hp as a health parameter in the milk or in the plasma / serum of the animals makes it possible to assess the general state of health. The proof could be implemented as an instrument of an animal welfare-oriented inventory management in the operational diagnostics. In the long run, this should make animal health management more organized and easier to carry out. Health data collected electronically in herd management programs and data and information platforms should be quickly and easily accessible. This could open up new possibilities for the objective comparison of farm animal health. For the Hp determination, there are several possibilities for sampling scenarios. The sampling and measurement can take place on an occasion-related basis or as a regular check.

Occasional measurements take place once or over a short period of a few days. The health markers should be recorded during differentiated, operational situations. Depending on the occasion of the measurement, different conclusions can be drawn from the results.

- Screening of the stock: Sampling of all animals once a month to determine overall health and well-being for a stock
  - serves as a basis for farm comparisons: stocks could be assessed due to variations in the baseline concentration of parameters and the proportion of animals outside the normal range.
  - Measurement in milk yield test samples is recommended.
- Screening of individual animals: occasion-related sampling of selected individual animals to determine the current state of health
  - Purchase / sale of cows: If an animal should leave the farm or be taken into the business, a health check can be carried out quickly and purposefully by measuring the health parameter in the milk.
  - Drying: Only udder-healthy cows may be kept dry. Undiscovered clinical mastitis or subclinical mastitis should be excluded.

For detected diseases, the measurements outside the milking routine could be continued on an ad hoc basis to monitor the therapy. In support of the puerperal control, the sensitive phase of early lactation could be monitored in the first few weeks, as the cows are particularly susceptible to disease due to the metabolic changes at the onset of lactation and the burden of calving.

### 2. INTENDED USE

---

The NovaTec VetLine Bovine Haptoglobin ELISA is intended for the quantitative determination of Haptoglobin in bovine serum or milk.

### 3. PRINCIPLE OF THE ASSAY

---

The quantitative immunoenzymatic determination of bovine Haptoglobin is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antibodies to bind corresponding Haptoglobin of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured Haptoglobin. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of Haptoglobin in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA microwell plate reader.

### 4. MATERIALS

---

#### 4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with anti-bovine Haptoglobin antibodies; in resealable aluminium foil.
- **Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015 % (v/v) CMIT/MIT (3:1). Can also be used as a negative control.
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 15 mL of peroxidase labelled anti-bovine Haptoglobin antibodies; coloured red; ready to use; white cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.

- **Standards:** 6 vials, each containing 2 mL standard (bovine Haptoglobin); coloured yellow; ready to use;  $\leq 0.02\%$  (v/v) MIT.
 

Standard A:	500	ng/mL; yellow cap
Standard B:	250	ng/mL; yellow cap
Standard C:	125	ng/mL; yellow cap
Standard D:	62.5	ng/mL; yellow cap
Standard E:	31.25	ng/mL; yellow cap
Standard F:	15.63	ng/mL; yellow cap

For hazard and precautionary statements see 12.1

For potential hazardous substances please check the safety data sheet.

## 4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

## 4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000  $\mu\text{L}$
- Vortex tube mixer
- Distilled water
- Disposable tubes

## 5. STABILITY AND STORAGE

---

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

## 6. REAGENT PREPARATION

---

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

### 6.1. Microtiterplate

The break-apart snap-off strips are coated with anti-bovine Haptoglobin antibodies. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

### 6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

### 6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

## 7. SAMPLE COLLECTION AND PREPARATION

---

Use bovine serum or milk samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

### 7.1. Sample Dilution

Before assaying, all samples should be diluted with Sample Dilution Buffer. Different dilutions have to be tested to obtain absorbance values within the calibration curve.

For serum samples we recommend to start with the following dilution: 1+50, 1+100, 1+1000, 1+10,000.

For milk samples we recommend to start with the following dilutions: 1+5, 1+10, 1+50, 1+100

For example: Dispense 10  $\mu\text{L}$  sample and 1 mL Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

## 8. ASSAY PROCEDURE

---

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to  $37 \pm 1$  °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour  $\pm$  5 min at  $37 \pm 1$  °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be  $> 5$  sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!  
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

### 8.1. Measurement

Adjust the ELISA microwell plate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA microwell plate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

**Measure the absorbance** of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

## 9. RESULTS

---

### 9.1. Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate Blank:** Absorbance value  $< 0.100$
- **Standard A:** Absorbance value  $> 0.800$
- **Standard B:** Absorbance value  $> \text{Standard C}$
- **Standard C:** Absorbance value  $> \text{Standard D}$
- **Standard D:** Absorbance value  $> \text{Standard E}$
- **Standard E:** Absorbance value  $> \text{Standard F}$
- **Standard F:** Absorbance value  $< 0.200$

**Standard A  $>$  Standard B  $>$  Standard C  $>$  Standard D  $>$  Standard E  $>$  Standard F**

If these criteria are not met, the test is not valid and must be repeated.

### 9.2. Calculation of Results

In order to obtain **quantitative results** plot the (mean) absorbance values of the 6 Standards A - F in a system of coordinates against their corresponding concentrations, starting with the lowest concentration (Standard F) (15,63 / 31,25 / 62,5 / 125 / 250 and 500 ng/mL) and draw a standard curve (absorbance values on the y-axis, concentrations on the x-axis).


Read results from this standard curve employing the (mean) absorbance values of each bovine sample.

For the calculation of the standard-curve mathematical Point to Point function should be used.

In order to calculate the final concentration, the values from the calibration line have to be corrected by the sample dilution.

**Example:** The sample was diluted 1+100 in Sample Dilution Buffer. The sample resulted in an OD of 0.579. According to the calibration line the concentration is e.g. 69.44 ng/mL. Now you have to multiply this with your dilution factor. In this example the factor is 100. For this reason the final concentration in the sample is 6944 ng/mL respectively **6.94 µg/mL**.

### 9.3. Typical Standard Curve

### 9.4. Interpretation of Results and Recommendations - Preliminary-

Normal value ranges for this ELISA should be established by each laboratory based on its own sample populations in the geographical areas serviced.

The following values should be considered as a guideline:

Bovine milk samples - Preliminary-		
Reference range [µg/mL]	Rating	Clinical status
< 0.5	low	healthy
1.0-4.5	critical	notable value
10-50	medium	subclinical mastitis
> 300	high	clinical mastitis
Bovine serum samples - Preliminary-		
10-50	low	healthy
600-700	medium	Physiological (5.-7.dpp*)
>1,000-2,000	high	sick
Diagnosis of a disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. * dpp- day post partum		

## 10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

### 10.1. Precision

-to be determined-

### 10.2. Diagnostic Specificity

-to be determined-

### 10.3. Diagnostic Sensitivity

-to be determined-

### 10.4. Analytical Sensitivity

The analytical sensitivity is defined as the apparent concentration of the analyte that can be distinguished from the zero calibrator.

-to be determined-

### 10.5. Interferences

-to be determined-

### 10.6. Cross Reactivity

-to be determined-

### 10.7. Measurement range

-to be determined-

## 11. LIMITATIONS OF THE PROCEDURE

---

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

## 12. PRECAUTIONS AND WARNINGS

---

- **For research use only!**
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

### 12.1. Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

#### Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray
P280	Wear protective gloves/protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet.

### 12.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Alsemgeest, S.P.M., Kalsbeek, H.C., Wensing, T., Koeman, J.P., van Ederen, A.M., Gruys, E., 1994. Concentrations of serum Amyloid-A (SAA) and haptoglobin (HP) as parameters of inflammatory diseases in cattle. *Veterinary Quarterly* 16, 21-23.

Eckersall, P.D., Connor, J.G., 1990. Plasma haptoglobin in cattle (*Bos taurus*) exists as polymers in association with albumin. *Comp. Biochem. Physiol. B* 96, 309-314.

Haghkhah, M., Nazifi, S., Jahromi, A., 2010. Evaluation of milk haptoglobin and amyloid A in high producing dairy cattle with clinical and subclinical mastitis in Shiraz. *Comparative Clinical Pathology* 19, 547-552.

Horadagoda, N.U., Knox, K.M., Gibbs, H.A., Reid, S.W., Horadagoda, A., Edwards, S.E., Eckersall, P.D., 1999. Acute phase proteins in cattle: discrimination between acute and chronic inflammation. *Vet. Rec.* 144, 437-441.

Huzzey, J.M., Duffield, T.F., LeBlanc, S.J., Veira, D.M., Weary, D.M., von Keyserlingk, M.A., 2009. Short communication: Haptoglobin as an early indicator of metritis. *Journal of Dairy Science* 92, 621-625.

Nazifi, S., Rezakhani, A., Koohimoghadam, M., Ansari-Lari, M., Esmailnezhad, Z., 2008. Evaluation of serum haptoglobin in clinically healthy cattle and cattle with inflammatory diseases in Shiraz, a tropical area in sothern Iran. *Bulgarian Journal of Veterinary Medicine* 11, 95-101.

Petersen, H.H., Nielsen, J.P., Heegaard, P.M.H., 2004. Application of acute phase protein measurements in veterinary clinical chemistry. *Veterinary Research* 35, 163-187.

Wenz, J.R., Fox, L.K., Muller, F.J., Rinaldi, M., Zeng, R., Bannerman, D.D., 2010. Factors associated with concentrations of select cytokine and acute phase proteins in dairy cows with naturally occurring clinical mastitis. *J. Dairy Sci.* 93, 2458-2470.

## **ABBREVIATIONS / ABKÜRZUNGEN / ABRÉVIATIONS / ABBREVIAZIONI / ABREVIACIONES / ABREVIATURAS**

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VetLine Bovine Haptoglobin ELISA

## Test Preparation

Prepare reagents and samples as described.  
 Establish the distribution and identification plan for all samples and standards/controls on the plate layout supplied in the kit.  
 Select the required number of microtiter strips or wells and insert them into the holder.

## Assay Procedure

	Substrate Blank (A1)	Standard A - F	Sample (dilution: refer to point 7.1)
Standard A - F	-	100 µL	-
Sample (dilution: refer to point 7.1)	-	-	100 µL
Cover wells with foil supplied in the kit <b>Incubate for 1 h at 37±1°C</b> Wash each well three times with 300 µL of Washing Buffer			
Conjugate	-	100 µL	100 µL
<b>Incubate for 30 min at room temperature (20...25 °C)</b> Do not expose to direct sunlight Wash each well three times with 300 µL of Washing Buffer			
TMB Substrate Solution	100 µL	100 µL	100 µL
<b>Incubate for exactly 15 min at room temperature (20...25 °C) in the dark</b>			
Stop Solution	100 µL	100 µL	100 µL
Photometric measurement at 450 nm (reference wavelength: 620 nm)			

## NovaTec Immundiagnostica GmbH

Waldstraße 23 A6  
 63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760

Fax: +49 (0) 6074-487629

Email: info@NovaTec-ID.com

Internet: www.NovaTec-ID.com

BHAPV4040-engl,dt-29042020-TL-RUO