# INGEZIM HAEMOPHILUS

Prod Ref: 11.HPS.K1

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## I. TECHNICAL BASIS

This is an immunoenzymatic assay for the detection of antibodies against Haemophilus parasuis (HPS) in porcine serum.

HPS is responsible for Glasse syndrome, a condition characterized by septicaemia which results in arthritis and meningitis. To date, 15 different serotypes have been identified with serotypes 1,5,10,12,13 and 14 being the most pathogenic.

The porcine serum samples and the controls are incubated in wells coated with HPS antigens (Ags) (odd columns), and in wells with a cell lysate (even columns) that serve as negative control. The antibodies (Abs) specific to HPS present in positive serum samples will bind to the Ags in the wells. After several washes to eliminate unbound

substances, a conjugate (an Ab coupled to an enzyme) targeted at the porcine Abs is added. After incubation, the excess of this conjugate is eliminated by a second wash and its attachment is revealed with a chromogenous substrate. Following this incubation, the enzyme, if present, reacts with the substrate and a green color develops. The reaction is stopped and the optical densities are read. The intensity of the colour allows the determination of the type of sample tested. A negative sample will show a weak reaction (pale green) whereas a strong positive will show a strong reaction (dark green). All shades of green between dark and pale represent various degrees of positivity.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
10. Substrate and stop solution must be handled with care.

## III. STORAGE OF COMPONENTS

All reagents and components must be stored between +2°C and +8°C. WARNING: Stop Solution may precipitate at this temperature, if it happens, warm up at +37°C before use.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLE

Samples and controls need to be tested as duplicate.

Dilute porcine serum samples at 1/50 in 1X serum diluent (e.g., 10µl sample in 490 µl 1x serum

diluent). Make sure you use a new tip for each sample. Also make sure each dilution is properly mixed before being distributed into the wells.

## VI. PREPARATION OF REAGENTS

- **Washing solution:**  
Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water (i.e. 40 ml of concentrate with 960 ml of water).
- **Preparation of Diluent:**  
Supplied diluent is 10x concentrated. Before using it, a 1/10 dilution must be done with distilled or deionized water (1 vol. of concentrated supplied Diluent with 9 volume of distilled water). When ready this solution remains stable between +2°C and +8°C.
- **Control sera:**  
Controls are supplied ready to use. Do not dilute
- **Preparation of the conjugate: to make immediately before use.**  
Dilute the needed quantity of conjugate provided in the kit 1/100 into diluent:  
→ The necessary and sufficient quantity for a complete plate is 110 µl of conjugate with 10,9 ml of diluent.  
Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.
- **Preparation of the substrate solution.**  
Substrate is supplied ready to use.

## VII. TEST PROCEDURE

Bring all reagents to room temperature and mix well manually before use.

1. Put strips coated with HPS Ags on the odds columns of a frame, and strips without Ags (blue) on the even columns. Put the rest of the strips in the plastic bags included for that purpose.
2. Make a schematic representation of the plate and the distribution of controls and samples.
3. Dispenses 100 µl ready-to-use positive control into wells A1 and A2.
4. Dispenses 100 µl ready-to-use negative control into wells B1 and B2.
5. Dispenses 100 µl diluted samples into wells C1/C2, D1/D2.
6. Incubate at  $23 \pm 2^\circ\text{C}$  for 30 minutes.
7. Wash each well with 4 times 300 µl 1X wash solution. (see section A). Throw away all liquid contained in the plate after each wash. After the last wash, dry the plate by tapping it on absorbent paper.
8. Dispense 100 µl diluted conjugate into each well.
9. Incubate at  $23 \pm 2^\circ\text{C}$  for 30 minutes.
10. Repeat step 7.
11. Dispense 100 µl ready-to-use substrate into each well.
12. Incubate, away from light, at  $23 \pm 2^\circ\text{C}$  for 20 minutes.
13. Dispense 100 µl ready-to-use stop solution into each well.
14. Measure optical densities (OD) at 405 nm. within 5 min after the addition of stop solution. If the microplate reader is equipped with a reference filter, set it a 490nm. The reading should be done no later than 15 minutes after the addition of the stop solution. Calculate the results.

## VIII. READING AND RESULT INTERPRETATION

For each sample (S) and control, subtract the OD obtained in the well without Ag (2) from the well containing Ag (1). To obtain **Ratio**, divide each result by the one obtained for the positive control (P)

$$\text{RATIO} = \frac{(\text{ODs1(Ag)} - \text{ODs2 (without Ag)})}{(\text{ODp1(Ag)} - \text{ODp2 (without Ag)})}$$

### Validation

The following criteria must be met with in order to validate the analysis:

- Negative control ratio must be less than 0.25
- Positive control corrected OD must be greater than 0.8

### Interpretation

Sample ratio less than 0.4 is considered negative  
Sample ratio greater or equal to 0.6 is considered positive.

Sample ratio less than 0.6 but greater or equal to 0.4 is considered suspicious.

It is recommended to collect a 2<sup>nd</sup> serum sample 2 weeks later and test it.

Developed and manufactured by:

Packaged and distributed in Spain by:  
INMUNOLOGIA Y GENETICA APLICADA, S.A.  
Av.de la Institución Libre de Enseñanza 39, 8ª  
planta  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

IT-73840  
IT-73780

ISO 14001 ISO 9001:2015  
9191.JING 9175.JING