# INGEZIM MOQUILLO IgM

Prod Ref: 15.CDM.K2

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## I. TECHNICAL BASIS

---

This kit is based on a capture enzymatic immunoassay. We make a brief description of the technique below:

A Monoclonal antibody (Mab) specific for dog IgM is fixed on a solid support (polystyrene plate). When a dog serum is added, the IgMs present in the sample are captured by the Mab adsorbed on the plate. After washing to eliminate all non fixed material from the serum sample, we add the viral antigen (recombinant nucleocapsid) which is captured by the IgM only if the dog serum sample contained specific immunoglobulin M to CDV.

After additional washing we add a peroxidase-conjugated Mab specific to CDV nucleocapside. After washing and addition of the substrate a colorimetric reaction will be developed which could be measured by an ELISA reader.

In this way the presence of colour means the presence of antibodies IgM against the virus in the dog sera, and the absence of colour the absence of specific antibodies.

The aim of this kit is to provide the users with a reliable diagnostic technique for this disease.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

---

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **IMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. Use a new tip for each serum sample.
6. For each utilisation of the kit, control positive and negative sera must be tested in a systematic way.
7. Stop solution is a strong acid. Handle with care.
8. Substrate is very sensitive to light and contamination.

## III. STORAGE OF COMPONENTS

---

All reagents and components must be stored at 4 °C.

## IV. INFORMATION ABOUT THE WASHING STEPS

---

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

Dilute the sera samples at **1/100** in kit diluent to be tested.

## VI. PREPARATION OF REAGENTS

### Washing solution :

Dilute one part of the concentrate washing solution provided in the kit into 9 parts of distilled or deionized water. When ready this solution remains stable when is stored between +2°C y +8°C.

### Positive and Negative Control:

Controls are ready to use. Do not dilute.

### Recombinant viral antigen:

Antigen is ready to use. Do not dilute.

### Preparation of the conjugate: to make immediately before use.

Dilute 1/100 with diluent.

- The necessary and sufficient quantity of conjugate for a complete plate is: 110 µl of conjugate in 11 ml of diluent.
- The necessary and sufficient quantity of conjugate for a 8 wells strip is: 10µl of conjugate in 1 ml of diluent.

Prepare only the quantity needed for each time because the remainder volume has to be rejected.

Shake very well the solution prior to use.

## VII. TEST PROCEDURE:

1. All reagents and wells to be used, must be allowed to come to room temperature before use.
2. Add 100 µl of positive and negative controls to two different wells. Add 100 µl of each dilution of sera samples to test, prepared following the previous instructions. We recommend running sample and controls, in duplicates. Seal the plate and **incubate for 15 minutes at 37°C**.
3. Wash 4 times following the procedure previously described.
4. Add 100 µl of antigen to each well. Seal the plate and **incubate 15 minutes at 37°C**
5. Wash 4 times following the described procedure
6. Add 100 µl of conjugate, prepared following the previous instructions, to each well. Seal the plate and **incubate for 15 minutes at 37°C**.
7. Wash 4 times following the described procedure.
8. Add 100 µl of substrate solution, to each well. Keep the plate for **5 min at room temperature**.
9. Add 100 µl of stop solution to each well.
10. Read the absorbances with an ELISA reader at 450 nm. Within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading must be done with an ELISA reader at **450 nm**.

### Kit Validation:

The test is valid when:

**O.D positive control > 0.8**

**O.D negative control < 0,15**

### Results Interpretation:

Using the positive control value, the cut off points values must be calculated as follows:

**Cut off = O.D 450 nm positive control x 0.2**

### NEGATIVE SAMPLES:

Samples with an OD lower than the Cut off value.

### POSITIVE SAMPLES:

In spite of the good sensitivity and specificity of this test, all the positive samples should be confirmed with other diagnostic method. Definitive diagnostic must not be based on a unique test result. These results must be correlated with other data: animal history, vaccination, and clinical symptoms.

Samples with an OD higher than the Cut off value. In this case you could distinguish two kinds of samples:

- low-medium titres (corresponding to IFI values of 1/20-1/80). These samples show OD values between positive cut off and 70% of OD positive control.
- high titres (corresponding to IFI values  $\geq$  1/160). These samples show OD values higher than 70% OD of positive control (0.7 x OD of positive control)

Developed and manufactured in Spain by:

EUROFINS INGENASA, S.A.  
C. de los Hermanos García Noblejas 39, 8ª  
planta  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: ingenasa@ingenasa.com  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in

by:

A large empty rectangular box with a blue border, intended for the distributor's name and contact information.