INGEZIM GLUTEN Quick
30.GL2.K2
COMPOSICION DEL KIT
KIT COMPOSITION

Reactivo Reagent	1 placa (1x8x12 pocillos) 1 plate box (1x8x12 wells)	
	Uni.	Vol.
Placas de microtitulación (8x12 pocillos) recubierta con un ACM de gliadinas, secalinas y hordeinas 96 well microtitration plates (12x8 wells) coated with a gliadin,secalin, and hordein-specific Mab	1	
Estandar europeo de gliadina a las concentraciones de 50 ng/ml, 25 ng/ml, 12,5 ng/ml, 6,25 ng/ml y 3,12 ng/ml listos para su uso Vials with Gliadin European Standard containing solutions of 50 ng/ml, 25 ng/ml, 12.5 ng/ml, 6.25 ng/ml and 3.12 ng/ml ready to use	5	1,5 ml
Control negativo, listo para su uso. Negative control, ready to use.	1	1.5 ml
AcM anti gliadina, secalina, hordeina, conjugado con peroxidasa, listo para su uso. Vials with an anti-gliadin, secalin and hordein Mab-peroxidase conjugate, ready to use.	1	15 ml
Tampón de extracción, listo para su uso Bottles with extraction buffer, ready to use.	1	110 ml
Solución de lavado (25 x concentrado) Bottles with Washing solution (25x concentrated)	1	65 ml
Diluyente (DE12-05) (5x concentrado) Bottles with Diluent (DE12-05) (5x concentrated)	1	65 ml
Sustrato (TMB) Bottles with Substrate (TMB)	1	15 ml
Solución de frenado (H ₂ SO ₄ 0,5M) Bottles with Stop Solution (H ₂ SO ₄ 0,5M)	1	15 ml

OTROS MATERIALES Y REACTIVOS NO INCORPORADOS EN EL KIT
 OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT:

Etanol
 Agua Destilada

Ethanol
 Distilled water

I. FUNDAMENTO TECNICO DEL KIT

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I. TECHNICAL BASIS

This kit is based on a double antibody sandwich enzyme immunoassay (DAS or capture Elisa). The technique is briefly described below:

The plate is coated with the **R5 monoclonal antibody (Mab)**, which specifically recognises a common epitope to gliadins, secalins and hordeins. Therefore, when a sample is added into the wells, the capture Mabs will bind to the existing gliadins in the sample. After a washing step a second R5 Mab specific to gliadin conjugated with Peroxidase is added.

This conjugated Mab will bind to the captured gliadins and after a further washing step the

peroxidase action will be detected by adding the TMB substrate. The change from a colourless substrate solution into a blue product and, after addition of stop solution, into a yellow colour, can be measured with an ELISA reader.

The content of gliadin in the samples will be determined by an interpolation of their OD values in the calibration curve, done with the gliadin standards supplied with the kit. The sensitivity of the assay is 3 ppm of gluten.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (22°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the kit reagents.
5. Do not use components after expiry dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
8. Use a new tip for each new sample.
9. Systematically include the controls and the standard curve each time the assay is run.
10. **IMPORTANT** The TMB Substrate must be handled with care as it is very sensitive to light and contamination: pipette from the substrate storage bottle a sufficient amount for the assay into a separate dark container prior to the colour development step.
11. The stop solution is a strong acid and therefore must be handled with care.

III. STORAGE OF THE KIT COMPONENTS:

- ✓ All components must be stored between +2°C and +8°C. **DO NOT FREEZE ANY OF THE KIT COMPONENTS.**
- ✓ Avoid exposure of the kit and the components to direct sunlight at any time, as some reagents are light sensitive.

IV. INFORMATION ABOUT THE WASHING STEPS:

The steps can be done using a squeeze bottle, an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 µl on each well. After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over, in order to avoid the possible exchange of contents from one well to another.
- Dispense a volume of 300 µl of washing solution in each well.
- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as are indicated in the kit instructions.
- Prior to emptying the content of the last washing step, verify that the next reagent to be added to the plate is ready to be used. Do not let the plate dry longer than strictly necessary.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

V. PREPARATION OF REAGENTS:

Washing solution:

Dilute one part of the concentrated washing solution provided in the kit with 24 parts of distilled or deionized water. Once prepared, this solution remains stable when stored at +4°C.

Dilution buffer:

Dilute one part of the concentrated dilution buffer provided in the kit with 4 parts of

distilled or deionized water. Once prepared, this solution remains stable when stored between +2°C and +8°C.

Conjugate, Negative Control and Standards:

They are supplied ready to use. Do not dilute. Add 100 µl of each reagent directly to the well.

VI. SAMPLES PREPARATION:

Extraction procedure for Heat-processed products or products containing Soya.

1. Weigh 0.25 g of a milled food sample and place it in a 10 ml tube.
2. Add 2.5 ml of the Extraction buffer to the tube containing the sample.
3. Close the tube tightly and cover the cap with parafilm to avoid evaporation.
4. Mix thoroughly by vortexing (5-10 seconds) and place the tubes in a rack.
5. Incubate the tubes between 20 and 60 minutes at room temperature. During this incubation period, vortex the sample 2-3 times.
6. Open the tubes, add 7.5 ml of 80% ethanol. Thoroughly disperse the samples by vortexing 10 to 60 seconds and incubate 20-60 minutes at room temperature **with agitation** (a rotatory shaker at 45 turns per minute is recommended).
7. Centrifuge the tubes for 10 min at 2000-2500g at room temperature
8. Transfer the supernatants to clean 10 ml tubes. The solution is now ready for running an ELISA.

Extraction procedure for Non Heat-processed products and non-containing Soya products.

1. Weigh 0.25 g of a milled food sample and place it in a 10 ml tube.
2. Add 10 ml of **60%** ethanol and incubate for 20-60 minutes at room temperature **with agitation**.
3. Proceed to step 7 of the *heat processed products extraction procedure*.

Extraction of samples containing tannins or polyphenols (chocolate, red wine, fruit juices...)

1. In a 10 mL propylene tube, weight 0.25 g of Gelatin and 0.1 g of Polyvinyl-pyrrolidone. (Polyvinylpyrrolidone Sigma PVP-360) (Fish Gelatine Sigma No. G-7041).
2. Weight 0.25 g of the milled food sample and add it to the tube.
3. Add 2.5 mL of Extraction Buffer
4. Close the tubes and seal them with parafilm in order to avoid evaporation due to heat.
5. Mix by vortex until the sample is completely mixed.
6. Incubate tubes between 20 and 60 minutes at room temperature.
7. Add 7.5 mL of 80% Ethanol to the sample and incubate between 20 and 60 minutes at room temperature **with agitation** (a rotatory shaker at 45 turns per minute is recommended).
8. Centrifugate the samples 10 minutes at 2000-2500 g at Room Temperature.
9. Collect the supernatant into clean 10 mL propylene tubes.
10. These extracts should be stored at Room Temperature.

EXTRACTION OF SWAB SAMPLES

IMPORTANT: *this procedure can only be used for a qualitative detection of gluten (absence/presence) in a surface sample.*

Swab samples should not be used for quantification, as the initial amount of sample in the swab is unknown.

Considering that using a swab only a little amount of sample can be taken, the extraction procedure is similar than the "normal" one, but using less volume of Extraction Buffer and Ethanol (as long as the original proportion is maintained):

1. Once the sample has been taken, cut the tip of the swab and put it into a 10mL propylene tube.
2. Add 1mL of Extraction Buffer to the tube.
3. Close the tube tightly with the cap and seal it using Parafilm in order to avoid evaporation.
4. Vortex for 5-10 seconds.
5. Incubate between 20 and 60 minutes at room temperature.
6. Add 3mL of Ethanol (80%), vortex for 10-60 seconds and between 20 and 60 minutes at room temperature **with agitation**. IMPORTANT: do not remove the swab from the tube at this step.
7. Remove the swab after the incubation and use the solution for the ELISA (dilutions 1/12.5 and 1/25).
8. Compare the OD values of both dilutions with the OD of the lower standard, in order to have an idea about the absence/presence of Gluten on the sample surface. The results should be

reported considering the Limit of Detection of the assay (3ppm of Gluten)

GENERAL INDICATIONS

The extracted samples may be kept 7 days at room temperature (or 1 month at 4°C), before performing the ELISA test. While stored, close tightly the vials with parafilm in order to avoid evaporation.

The samples should be diluted as follows before adding to the ELISA plate wells:

- 1) Samples with estimated gluten contents lower than 10 ppm: make a 1/12.5 dilution, using the supplied diluent.
- 2) Samples with gluten contents between 10-100 ppm: make 1/25 dilutions.
- 3) Samples with gluten contents over 100 ppm: make 1/50 dilutions.

Dilution procedure:

1/12.5 = 920µl of diluent+80µl of sample

1/25 = 960µl of diluent+40µl of sample

1/50 = 980µl of diluent+20µl of sample

In general, for samples with unknown gluten content, we recommend a 1:25 dilution

VII. TEST PROCEDURE:

1. Prior to starting the test, equilibrate all reagents to room temperature (22-25°C).
2. Add 100 µl of sample dilutions. We recommend the use of two wells per diluted sample. Add 100 µl of standards and negative control. Seal the plate **and incubate 20 minutes at room temperature (22-25°C)**.
3. Wash 3 times following the procedure previously described in point IV (Washing steps).
4. Using a multi-channel pipette, add 100 µl of the ready to use conjugate to each well. Seal the plate **and incubate 20 minutes at room temperature (22-25°C)**.
5. Wash 3 times following the procedure previously described.
6. Using a multi-channel pipette, add 100 µl of TMB substrate to each well. Keep the plate for **20 min at room temperature**
7. Add 100 µl of the stop solution to each well in the same order in which the substrate was added. A positive reaction will change the colour from blue to yellow.
8. Read the OD of each well at 450 nm within the following 5 min after the addition of the stop solution.

VIII. READING AND INTERPRETATION OF THE RESULTS:

Validation of the test

The test can be considered valid when:

- ⇒ The OD of the 50 ng/ml standard point is higher than 1.5 and the OD of the Negative Control is lower than the 3,12 ng/ml point.

Interpretation:

The OD of the sample must be calculated as the arithmetic mean of the duplicated OD values.

Plot the calibration curve using the mean of the OD values obtained for the five standards (abscises) and their corresponding ng/ml values (ordinate).

Determine the concentration of gliadin in the sample (C) by interpolating the absorbance value of the sample in the curve.

The content (ppm) of gluten in a sample will be calculated as follows:

$$\text{ppm of gluten} = C \times D \times 2 \times 40 / 1000$$

where:


- C:** is the concentration of gliadin in the sample calculated from the calibration curve in ng/ml.
- D:** is a dilution factor of the sample (12.5, 25, 50)
- 2:** is a factor applied in order to express the results as concentration of gluten (instead of gliadin)
- 40:** is the dilution factor applied in the sample preparation (0,25g in 10 ml of extraction solution).
- 1/1000:** is the conversion rate from ng/ml to ppm

ADDITIONAL INFORMATION

- To ensure that the calculations obtained are correct, the OD values of the samples should be within the range of the OD values of the standard curve. It may be necessary to repeat the test using different dilutions of the samples to achieve these results. If so, it is important to apply the dilution factor in the calculation of results order to obtain correct values.
- Only in the case that the OD value of the sample is slightly higher than the OD of the 50 ng/ml point (up to a 15%), the assay can be considered as valid.

EXAMPLE:


- 1) OD of Negative Control = 0,13
- 2) OD of the standards points:

ng/ml	Abs450
50	2.768
25	1.660
12.5	1.039
6.25	0.717
3.12	0.361

- 3) The concentration of gliadin (C) in the sample can be obtained by extrapolating the OD value of the sample in the standard curve or by using an appropriate software. For example:

- ✓ Software: Excel (Microsoft)
- ✓ Once obtained the gliadin standards OD values, click on "Graphic assistant"
- ✓ Select type of chart: "XY (dispersion)"
- ✓ Click on the top window "chart" and then select "Add tendency line"
- ✓ In the tab "type", click on "polynomial"
- ✓ In the tab "Options", click on "Display the formula in the chart"
- ✓ In this equation, replace the unknown factor "x" by the sample OD value.

Sample	Dilution(D)	DO	ng/ml (C)*	Ppm gluten
1	12.5	0.21	< 3.12	< 3
2	12.5	0.500	4	4
3	25	0.893	10	20
4	25	2.817	51.5	103
5	50	2.532	44	176

* Obtained by replacing "x" by the OD value in the tendency line.

OVERVIEW OF PROCEDURE

1. Reagents must be equilibrated at room temperature (**22-25°C**) prior to starting the test.
2. Add 100µl of the **samples** (prepared and diluted as described in point VI) in duplicate wells.
3. Add 100µl of **standards** and **negative control**. Seal the plate with adhesive foil.
4. Incubate at room temperature (**22-25°C**) for 20 min.
5. Wash the wells 3 times with the *diluted wash* solution to remove unbound reagents.
6. Add 100µl of *Mab Peroxidase Conjugated* to each well. Seal with adhesive foil.
7. Incubate at room temperature (**22-25°C**) for 20 min.
8. Wash the wells 3 times with the *diluted wash* solution to remove nonreactive conjugate.
9. Add 100µl of **substrate (TMB)** to each well.
10. Incubate **20min** at room temperature
11. Add 100µl of stop solution to each well.
12. Read at 450nm.

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INMUNOLOGIA Y GENETICA APLICADA, S.A.