INgezim® Neospora 3.0

R.12.NC.K1

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INgezim® Neospora 3.0

R.12.NC.K1

INgezim® Neospora 3.0 is an immunoenzymatic assay based on an indirect ELISA technique, which uses a monoclonal antibody (Mab) specific to bovine IgG and a soluble extract of *Neospora caninum* proteins.

KIT FEATURES

APPLICATION


Detection of specific antibodies to *Neospora caninum* in bovine sera samples.

TECHNICAL BASIS OF THE KIT

1. Plates supplied are coated with *Neospora caninum* antigen. Samples are added to the wells and incubated.
2. If the sample contains antibodies to *Neospora caninum*, they will bind to the antigen.
3. When a monoclonal antibody specific for bovine IgG is added, it will bind to them. The monoclonal antibody has been conjugated with HRPO for future detection.
4. This binding is detected by the development of the colorimetric reaction after the addition of the substrate.

INTERPRETATION OF THE RESULTS

The assay establishes two cut offs. Samples will be considered **Positive** if its IRPC value is higher than the positive cut off; **Negative**, if its IRPC value is lower than the negative cut off. IRPC values between both cut offs will be considered **Doubtful**.

ASSAY VALIDATION

1. COMPARISON WITH OTHER TECHNIQUES

A comparative study was carried out with the previous version in order to evaluate the improvement of version 3.0. For this purpose, 29 sera positive by ELISA "in house" and Western Blot (WB), which is considered in the study as the "Gold Standard"; 11 sera negative by WB and ELISA "in house" and 7 sera positive by WB and negative by ELISA "in house" were used.

The following results were obtained:

| | | ELISA "in house" | WB |
|----------------------|-------------|------------------|------|
| INgezim Neospora 3.0 | Sensitivity | 93% | 81% |
| | Specificity | 89% | 100% |

2. COMPARATIVE STUDY USING CIVTEST® BOVIS NEOSPORA

Sera previously classified

A comparative study was performed with the commercial ELISA assay CIVTEST® BOVIS NEOSPORA.

29 sera positive by WB and ELISA "in house" (group A), 11 sera negative by WB and ELISA "in house" (group B), 7 sera positive by WB and negative by ELISA "in house" (group C), 33 sera not infected with *N. caninum* and negative by ELISA (group D) and 11 sera from animals infected with other protozoa forming cysts. (Besnoitia and/or Sarcocystis) (group E) were used.

The following results were obtained:

| | | ELISA "in house" | WB | CIVTEST® Bovis Neospora |
|-------------------------|-------------|------------------|------|-------------------------|
| CIVTEST® Bovis Neospora | Sensitivity | 93% | 81% | |
| | Specificity | 97% | 100% | |
| INgezim® Neospora 3.0 | Sensitivity | 93% | 81% | 100% |
| | Specificity | 97% | 100% | 100% |

*Comparative with WB was made only with groups A, B, & C. Groups D, and E were used together with the rest groups only to compare the 3 ELISAs.

Field sera

In this study, 614 field sera from areas in which the disease is present, were used. All sera were analyzed by INgezim® Neospora 3.0 and CIVTEST® Bovis Neospora.

Results obtained are shown below:


| | | CIVTEST® BOVIS NEOSPORA | | |
|----------------------|----------|-------------------------|----------|-------|
| | | POSITIVE | NEGATIVE | TOTAL |
| INGEZIM NEOSPORA 3.0 | POSITIVE | 63 | 3 | 66 |
| | NEGATIVE | 10 | 538 | 548 |
| | TOTAL | 73 | 541 | 614 |

From this study, we can conclude that, in comparison with CIVTEST® Bovis Neospora, INgezim® Neospora 3.0 showed 86% and 99% relative sensitivity and specificity respectively.

3. ROC ANALYSIS

For this analysis, all sera classified by CIVTEST® Bovis Neospora used in the validation, were included (studies 1 and 2).

The obtained results indicated a good capability of INgezim® Neospora 3.0 to discriminate between positive and negative samples.

COMPOSITION OF THE KIT

- 96 wells microtitration plates
- Vials with positive control
- Vials with negative control
- Vials with conjugate
- Bottle with washing solution
- Bottle with diluent
- Bottle with substrate
- Bottle with stop solution

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Eurofins -INGENASA, S.A