INGEZIM INFLUENZA PORCINA 2.0

Prod Ref: 11.FLU.K1

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11.FLU.K1

I. TECHNICAL BASIS

This kit is based on an indirect enzyme immunoassay (Indirect ELISA). A brief description of the technique is given below:

The antigen is coated on polystyrene plates. When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on the plate. After washing to eliminate all non-fixed material the presence of swine antibodies can be detected using a

specific peroxidase conjugate. After addition of the substrate a colorimetric reaction will appear that can be measured by a spectrophotometer.

The presence of colour indicates the presence of antibodies against the virus in the swine sera, and the absence of colours, the absence of specific antibodies.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Substrate must be handled with care, it is very sensitive to light and contamination
10. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
11. Stop solution is a strong acid. Handle with care

III. STORAGE OF COMPONENTS

Plates and reagents must be stored between +2°C y +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF REAGENTS

• **Washing solution:**

Dilute one part of the concentrated washing solution provided in the kit with 24 parts of distilled or deionized water (40 ml of concentrated solution and 960 ml of water). Once ready this solution remains stable at +4°C.

• **Control Sera, Diluent and Conjugate:**

All these components are ready to use. Do not dilute.

VI. PREPARATION OF SAMPLES:

The sera samples need to be at 1/100 dilution (i.e. 5µl of serum with 495 µl of diluent)

VII. TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. Add samples diluted as it's specified in the previous instructions to the wells of the plate. Add 100 µl of positive and negative control serum. We recommend running samples and controls in duplicate for confirmatory purposes. Incubate the plate for 1 hour at 37°C ± 1°C
3. Wash 4 times following the described procedure.
4. Add 100 µl of conjugate to each well. Incubate the plate for 1 hour at room temperature (+18-25°C).
5. Wash 5 times following the described procedure.
6. Add 100 µl of substrate to each well. Keep the plate for 10 min at room temperature.
7. Add 100 µl of stop solution to each well.
8. Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at 450 nm.

A. Validation of the test:

The test is considered valid when:

- The ratio:

$$\frac{\text{OD Negative Control}}{\text{OD Positive Control}} < 0.2$$

- OD of positive control ≥ 1

B. Interpretation of the results:

Calculate the ratio: S/P

$$\frac{\text{Sample OD}}{\text{Positive control (OD)}}$$

- Samples with S/P values ≥ 0.2 are positive for antibodies to influenza A viruses.
- Samples with S/P values < 0,2 will be considered negatives in antibodies to Influenza A viruses

C. Titer of the samples:

The titer of a sample must be calculate as follow:

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S.A.