

## INgezim Brucella Ovis

R.13.B0.K1

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

## INgezim Brucella Ovis

R.13.BO.K1

**INgezim Brucella ovis** is an immunoenzymatic assay based on indirect ELISA technique, which uses a monoclonal antibody (MAb) specific for ruminant's IgG immunoglobulins, and as antigen, Hot Saline (HS) extract of *B. ovis*.

### TECHNICAL BASIS OF THE KIT

1. Plates are coated with *Brucella ovis* antigen (HS). Samples are added to the wells and incubated.
2. If the sample contains antibodies to *Brucella ovis*, they will bind to the antigen.
3. To detect this reaction a MAb-PO specific to small ruminant IgG is added. This binding is detected by the development of a colorimetric reaction after the addition of the substrate.

### APPLICATION

Detection of specific antibodies to *Brucella ovis* in ram serum samples.

### INTERPRETATION OF THE RESULTS

For the results interpretation, one cut off is used: Samples will be considered **Positive** if their OD values are equal to or higher than the cut off; **Negative**, if their OD values are lower than the cut off.

## VALIDATION

### Diagnostical Sensitivity and Specificity:

In order to evaluate the performance of the assay, 1499 ram's serum samples were analyzed:

- 338 *B. ovis* isolation negative.
- 203 *B. ovis* isolation positive.
- 704 from herds that were infected with low prevalence and then became negative and from herds historically negative.
- 230 Complement Fixation (CF) negative to *B. ovis*.
- 24 CF positive to *B. ovis*.

**Relative sensitivity respect to isolation 98,5%**

**Relative sensitivity respect to CF >99%**

**Relative specificity respect to free status 98%**

**Relative specificity respect to CF 98,7%**

When the data of these 1499 sera were analyzed as a total, **99% sensitivity and 98% specificity** were obtained.

	Positive to <i>B. ovis</i>	Negative to <i>B. ovis</i>	TOTAL
ELISA +	225	25	250
ELISA -	2	1247	1249
TOTAL	227	1272	1499

### COMPOSITION OF THE KIT

- Microtitration plates of 96 wells divided into 12 strips of 8 wells each
- Vial with Positive Control.
- Vial with Negative Control.
- Vial with Conjugate.
- Bottles with Washing Solution.
- Bottles with diluent.
- Bottle with substrate (TMB) ready to use.

PRODUCT MANUFACTURED BY INGENASA  
REGISTRATION NUMBER: 3375 RD

IT-73840 ISO 14001:2015 ISO 9001:2015  
IT-73780 9191.INGE 9175.ING2

Shelf life: **18 months**  
Store at 2°C-8°C

Ed.020217