
INgezim ARTERITIS 2.0

R.14.EA2.K1

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
INgezim ARTERITIS 2.0

R.14.EA2.K1

INgezim ARTERITIS 2.0 is an indirect immunoenzymatic assay which uses a monoclonal antibody specific of equine IgG and a negative and a positive antigen.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with negative and positive antigens of Equine Arteritis virus (EAV). Serum samples are added to the wells and incubated.
2. If the sample contains specific antibodies to EAV, these antibodies will bind to the positive antigen coating the plate.
3. Once the MAb specific of equine IgG conjugated with PO is added, it will bind to the antibodies present in serum, in the case of positive animals. This binding is detected by the development of a colorimetric reaction after the substrate addition.

APPLICATION

Detection of specific antibodies to Equine Arteritis virus.

INTERPRETATION OF THE RESULTS

The assay defines two cut offs. Samples are considered **positive**, if $OD Ag_{POS} - OD Ag_{NEG} > 0.15$. Samples are considered **negative**, if $OD Ag_{POS} - OD Ag_{NEG} < 0.1$

SENSITIVITY

The sensitivity of the assay was determined using 87 sera of animals, positive by Seroneutralization (SN). **The results obtained indicated 98.8% sensitivity respect of SN.**

SPECIFICITY

180 sera negative by SN were used to determine the specificity of the assay. **The results obtained indicated 96.6% specificity respect of SN.**

ANALYTICAL SENSITIVITY

The OIE Reference sera for SN were analyzed.

- Strong positive reference serum
- Medium positive reference serum
- Weak positive reference serum
- Negative reference serum

The results obtained indicated that the assay **fulfils the OIE requirement concerning sensitivity and specificity.**

COMPOSITION OF THE KIT

- Microtitration plates of 96 well coated with negative antigen
- Microtitration plates of 96 well coated with positive antigen
- Vials with Positive Control
- Vials with Negative Control
- Vials with Peroxidase Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with stop solution
- Bottle with substrate (TMB) ready to use

PRODUCT MANUFACTURED BY INGENASA
REGISTRATION 0606 RD

SHELF LIFE: 18 months
Stored at 2°C-8°C