

INgezim CORONA FELINO

R.16.FCV.K1

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

INgezim CORONA FELINO

R.16.FCV.K1

Ingezim CORONA FELINO is based on an indirect enzymatic immunoassay (indirect ELISA), which uses a specific monoclonal antibody (MAb), specific to feline immunoglobulins (Igs) and a Feline Infectious Peritonitis Virus (FIPV) antigen.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with antigen FIPV. Serum samples are added and incubated.
2. If the samples contain specific antibodies to either of the two feline Coronavirus, they will bind to the antigen.
3. When a MAb-PO specific of feline immunoglobulins is added, it will bind to the Igs previously bound to the antigen. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

- Detection and/or titration of specific antibodies to feline Coronavirus in cat serum samples.
- An approximated differential diagnosis between Feline Infectious Peritonitis (FIPV) and Feline Enteric Coronavirus (FECV) can be made based on the titre of each sample. Estadistically 75% of FIPV infected cats have titres equal to or higher than 1/6400 while only 20% of FECV affected cats show such high titres.

INTERPRETATION OF THE RESULTS


Screening: OD values higher than the cut off are considered **Positive**; OD values lower than the cut off are considered **Negative**.

Titration: The titre of each sample will be the last dilution at which the OD value is higher than the cut off.

DIFFERENTIAL DIAGNOSIS BASED ON TITRES

The study was made using 109 feline sera: 29 positive to FIPV and 80 to FECV previously classified by other assays. The results obtained are shown below:

FECV INFECTED			FIPV INFECTED		
TITER	Nº SERA	PERCENTAGE	TITER	Nº SERA	PERCENTAGE
100	4	5,0%	100	0	0%
200	5	6,2%	200	0	0%
400	7	8,7%	400	0	0%
800	11	13,7%	800	3	10%
1600	20	25,0%	1600	2	7%
3200	16	20,0%	3200	3	10%
6400	11	13,7%	6400	6	20%
12800	4	5,0%	12800	5	17%
25600	2	2,5%	25600	10	34%

COMPOSITION OF THE KIT

- 96 wells microtitration plates
- Vials with Positive Control
- Vials with Negative Control
- Vials with conjugate
- Bottle with washing solution
- Bottle with diluent
- Bottle with substrate
- Bottle with stop solution

PRODUCT MANUFACTURED BY INGENASA

Shelf life: **18 months**
Store at: 2°C-8°C

Ed. 020217