CPV

## vp2 CPV specific Monoclonal Antibody Anticuerpo Monoclonal específico de vp2 de CPV

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*Canine Parvovirus (CPV) belongs to the feline parvovirus subgroup of the genus Parvovirus, from the Parvoviridae family. The CPV capsid is icosahedral and is composed of 60 copies of structural proteins, mostly VP2 and some VP1. The VP2 is the N-terminal truncated form of VP1. When some VP2 molecules undergo an N-terminal proteolysis, they produce the VP3. This virus is one of the major new pathogens of dogs and has caused a worldwide pandemic. The infection with Canine Parvovirus may show very different clinical signs, ranging from generalized neonatal infections to enteritis, cerebral hypoplasias, myocarditis, etc. The virus is secreted from day 2 post infection today 8 p.i. For diagnosis confirmation, the most reliable way to confirm the disease is the virus isolation from faeces. Occasionally, the detection of specific antibodies to CPV can be a useful tool from a diagnosis and epidemiological point of view (optimal moment for vaccination, etc)*

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*Detection of CPV by immunofluorescence (IFA) and immunoblotting (IB)*

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
*The hybridome which produces the monoclonal antibody has been obtained by the fusion of lymphocytes from Balb/c mice's spleen with myeloma X63/Ag8653 cells. The IgG has been purified by affinity chromatography, showing a purity of 99%.*

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*Briefly, the cells were fixed with a 50% (v/v) mixture of methanol-acetone at -20°C. They were then rehydrated and incubated with the corresponding MAb (See recommended use conditions for dilution. Note that the supernatant is supplied ready to use) for 60 min. at room temperature, washed with Wash with PBS-0.25% BSA and then stain with fluorescein-labelled rabbit anti-mouse antibodies for another period of 45 min. Finally, slides were washed with PBS, dried and set with 90 % glycerol-PBS.*

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