

INgezim LEISHMA-CROM (Ab)

R.15.LSH.K41

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

INgezim LEISHMA-CROM (Ab)

R.15.LSH.K41

INgezim LEISHMA-CROM is based on an indirect immunocromatography assay technique, which uses a monoclonal antibody (MAb) specific of canine IgG, and an inactivated antigen of *L. infantum*.

TECHNICAL BASIS OF THE KIT

The device consists of a plastic casing with three windows:

- Window1: Where the Monoclonal Antibody specific of canine IgG and the Control Protein are conjugated to colloidal gold particles.
- Window2: Where the inactivated antigen of *L. infantum* is adsorbed.
- Window 3: Where the MAb specific to the control protein is attached to indicate the validation of the assay.

APPLICATION

Detection of antibodies specific of Leishmania in canine sera or plasma samples.

INTERPRETATION OF THE RESULTS

Positive sample: Red/pink lines in the Test and Control zones.
Negative sample: Only a red/pink line in the Control zone. **Not valid assay:** No line in the Control zone

VALIDATION

Correlation with ELISA

Study 1.

43 canine field sera (20 positive, 18 negative, 5 doubtful) previously classified by INgezim Leishmania, were analyzed. Considering doubtful ones as positive, the results indicated **95%** correlation between assays.

		ELISA (INGEZIM LEISHMANIA)		
		POSITIVE	NEGATIVE	TOTAL
INGEZIM LSHCROM	POSITIVE	23	0	23
	NEGATIVE	2	18	20
	TOTAL	25	18	43

Study 2.

86 canine field sera (45 positive, 41 negative) previously classified by INgezim Leishmania, were analyzed. The results indicated **96.5%** correlation between assays.

		ELISA (INGEZIM LEISHMANIA)		
		POSITIVE	NEGATIVE	TOTAL
INGEZIM LSHCROM	POSITIVE	42	0	42
	NEGATIVE	3	41	44
	TOTAL	45	41	86

COMPOSITION OF THE KIT

- Diagnostic devices individually wrapped
- Diluent tubes, to prepare the dilution of the sera samples
- Dropper bottle dilution buffer
- Disposable tips

PRODUCT MANUFACTURED BY INGENASA
REGISTERED 1110 RD

SHELF LIFE: **24 months**
Stored: 4°C-25°C

Ed. 020217