

INGEZIM INFLUENZA PORCINA

R.11.FLU.K1

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

INgezim INFLUENZA PORCINA

R.11.FLU.K1

INgezim Influenza Porcina is based on a Indirect ELISA technique which uses a monoclonal antibody (MAb) specific of porcine IgG.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with porcine Influenza antigen. Serum samples are added and incubated.
2. If the samples contain specific antibodies to Porcine Influenza Virus, they will bind to the antigen.
3. When the MAb-PO specific of porcine immunoglobulin is added, it will bind to the Igs bound to the antigen. The binding is detected by the development of a colorimetric reaction after the addition of the substrate.

APPLICATION

- o Detection and/or titration of specific antibodies to Influenza virus in porcine sera samples.
- o Vaccination control.
- o Valid for strains H₁N₁, H₁N₂ y H₃N₂.

INTERPRETATION OF THE RESULTS


The assay uses one Cut off. Samples showing OD values higher than or equal to the cut off will be considered **positive**. Samples showing OD values lower than the Cut off must be considered **negative**.

APPLICATION

Studies with two groups of 7 pigs were made. Pigs were vaccinated with different doses of a commercial vaccine which contains antigen of strains H₁N₁, H₁N₂ and H₃N₂. Extractions were taken at days 0, 28, 42 and 49 and were analyzed to determine the performance of the assay detecting specific antibodies to the three strains and the precocity of detection.


1. Detection of response to different strains

The results obtained indicated that the assay is able to detect antibodies specific of the three strains of Influenza affecting pigs (H₁N₁, H₃N₂, H₁N₂). The figure represents some examples.

2. Early detection of specific antibodies

The results obtained with different extractions made at different days post vaccination were analysed, showed that the assay is able to detect specific antibodies to Influenza at day 28 post vaccination. The behaviour is similar for the three strains affecting pigs.

3. CORRESPONDENCE WITH IHA

24 samples of animals inoculated with strains, H₁N₁, H₁N₂, H₃N₂ or a mix of these 3 ones were analyzed. These samples had been previously tested by IHA. Results obtained showed 87% sensitivity and 89% specificity, being the 2 false negative samples near the cut off and with low titers of IHA (1/20 & 1/40).

COMPOSITION OF THE KIT

- Microtitration plates of 96 wells
- Vials with Positive control
- Vials with Negative control
- Vials with Conjugate
- Bottle with washing solution
- Bottle with diluent
- Bottle with substrate
- Bottle with stop solution

PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: **15 months**
Store at 2°C-8°C

Ed 280917