INGEZIM GLUTEN SEMIQ

R.30.GSQ.K2

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
INGEZIM GLUTEN SEMIQ

R.30.GSQ.K2

Ingezim gluten is a semi quantitative immunoenzymatic assay which uses the R5 monoclonal antibody specific of Prolamins. **Limit of Detection: 10 ppm.**

TECHNICAL BASIS OF THE KIT

1. 4 strips plates are supplied coated with a **monoclonal antibody (MAB) specific for Gliadin, Secalin and Hordein**. Samples are added to the wells and incubated for 15 minutes.
2. If the sample contains the antigen, it will bind to the antibody specific for Gliadin coating the plates.
3. When MAb-PO specific of Gliadin is added, after 15 minutes incubation, it will bind to the antigen of the sample.
4. Binding is detected by a colorimetric reaction after the addition of the substrate during 5 minutes


INTERPRETATION OF THE RESULTS

Comparison between the OD value of the sample and the OD value of the SemiQ point

Assay endorsed by the CODEX Alimentarius as Type I Method for the determination of gluten in food samples

VALIDATION

Comparing with the reference kit INGEZIM GLUTEN

SENSITIVITY

The sensitivity of the assay is 10ppm of gluten (5ppm of gliadin)

SPECIFICITY

The specificity of the MAB used in the assay was determined by the immunoblotting technique. The MAB is able to detect specifically the prolamins of gluten present in WHEAT, RYE and BARLEY (gliadin, secalin and hordein respectively).

COMPOSITION OF THE KIT

- Microtitration plate of 32 well
- Vial with SemiQ Point
- Vial with PO-Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with substrate (TMB) ready to use
- Bottle with stop solution
- Bottle with Extraction Solution

PRODUCT MANUFACTURED BY INGENASA

SHELF LIFE: 10 months
Stored at 4°C

02/02/2018