
INGEZIM GLUTEN Quick

R.30.GL2.K2

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
INGEZIM GLUTEN Quick

R.30.GL2.K2

Ingezim gluten Quick is a RAPID immunoenzymatic assay which uses the R5 monoclonal antibody specific of Prolamins.
Limit of Detection: 3 ppm.

TECHNICAL BASIS OF THE KIT

1. Plates are supplied coated with **R5 monoclonal antibody (MAb) specific for Gliadin, Secalin and Hordein**. Samples are added to the wells and incubated.
2. If the sample contains the antigen, it will bind to the antibody specific for Gliadin coating the plates.
3. When MAb-PO specific of Gliadin is added, it will bind to the antigen of the sample.
4. Binding is detected by a colorimetric reaction after the addition of the substrate

APPLICATION

The test has been designed for samples with a content of Gluten between 0 – 300 ppm

The assay is designed for performing short incubations (20 minutes), for the extraction and analysis.

Assay endorsed by the CODEX Alimentarius as Type I Method for the determination of gluten in food samples

VALIDATION

REFERENCE SAMPLES

Different samples (Spike samples, FAPAS samples and samples with a known content of Gluten) were analysed by both kits: INGEZIM GLUTEN and INGEZIM GLUTEN QUICK.

All the samples were analysed three times by different people. The average results are showed in the table below

	Theoretical ppm	Gluten Quick	Gluten
Spike 0 ppm	0	<1,56	<1,56
Spike 20 ppm	20	19	24
Spike 50 ppm	50	50	58
Spike 100 ppm	100	104	114
Spike 200 ppm	200	216	196
Fapas Cake Mix	20.5	22	22
Fapas 20 ppm	20	22	22
M 104 65 ppm	65	66	54
10 ppm A+	10	13	13
100 ppm A+	100	106	105
M 83 21 ppm	21	22	29
M 93 25 ppm	25	26	28
Fapa 2795A	70	71	81

ANALYTICAL SENSITIVITY

The sensitivity of the assay is 3ppm of Gluten (1.5 ppm of Gliadin).

SPECIFICITY

The specificity of the MAb used in the assay was determined by the immunoblotting technique. The MAb is able to detect specifically the prolamins of gluten present in WHEAT, RYE and BARLEY (gliadin, secalin and hordein respectively).

KIT COMPOSITION

- Microtitration plate of 96 well
- Vial with Negative Control
- European Gliadin Standards (5)
- Vial with PO-Conjugate
- Bottle with Washing Solution
- Bottle with diluent
- Bottle with substrate (TMB)
- Bottle with stop solution
- Bottle with Extraction Solution

PRODUCT MANUFACTURED BY INGENASA

IT-73840 ISO 14001:2015 ISO 9001:2015
IT-73780 9191.INGE 9175.ING2

SELF LIFE: **10 months**
Stored at 4°C

02/02/2018