
INGEZIM OVOALBÚMINA

R.30.OVO.K2

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
INGEZIM OVOALBÚMINA

R.30.OVO.K2

Ingezim ovoalbumina is an immunoenzymatic assay which uses specific antibodies of ovalbumin (OVO), the main protein in the egg white. **Limit of Detection: 0.25 ppm.**

TECHNICAL BASIS OF THE KIT

1. Plates are supplied coated with Antibodies specific for OVO. Samples are added to the wells and incubated.
2. If the sample contains the antigen, it will bind to the antibody specific for OVO coating the plates.
3. When Ab-PO specific of OVO is added, it will bind to the antigen of the sample.
4. Binding is detected by a colorimetric reaction after the addition of the substrate

APLICACION

This assay is a rapid and reliable test that significantly reduces the time required to screen food products and environmental samples for the presence of ovalbumin residues.

Results can be read visually (qualitative interpretation), or using an ELISA reader (quantitative interpretation).

SENSITIVITY

The assay is able to detect **0.25ppm** of ovalbumin, by using a linear regression curve. It leads to more precise results (mainly at low concentrations).

Regression curve

RECOVERY RATES

1. IN RED WINE SAMPLES

Twelve red wine samples were spiked at 1ppm of ovalbumin. Performing 2 different extractions, the average recovery rate was **87%**.

2. IN SOLID FOODSTUFFS

Several food matrixes (pasta, rice flour, baby food, bread, chocolate cookies, custard...) were spiked at 1 ppm of ovalbumin. The average recovery rate was **104%**.

KIT COMPOSITION

- Microtitration plate of 48 wells
- Vial with Negative Control
- OVO Standards (4)
- Vial with PO-Conjugate
- Bottle with Washing Solution
- Bottle with substrate (TMB)
- Bottle with stop solution
- Bottle with Extraction Solution

PRODUCT MANUFACTURED BY INGENASA

SELF LIFE: **12 months**
Stored at 4°C

02/02/2018