INGEZIM BLV CONFIRMATION

Prod Ref: 12.BLV.K1

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12.BLV.K1

I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). The technique is briefly described below:

We fix positive and negative antigens on different plates. When a sample serum contains specific antibodies against the virus, they will bind to the viral antigen adsorbed on the plate. After a washing step to eliminate all non-fixed material from the sera sample, the presence of bovine immunoglobulins can be

detected using a specific peroxidase conjugated. After the addition of the substrate a colorimetric reaction will appear which can be measured by an ELISA reader.

This way the difference between the OD obtained with the positive and the negative antigen will determine the presence or absence of specific antibodies against the virus in the bovine serum.

The aim of this kit is to provide the users with a reliable and rapid confirmatory diagnostic technique for this disease.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the kit's reagents.
5. Do not use components after expiration dates and do not mix components from different batches.
6. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Systematically include a positive control and a negative serum each time the assay is run.
10. Substrate must be handled with care, as it is very sensible to light and contamination.
11. Stop solution is a strong acid solution and therefore must be handled with care. In case of accidental contact with skin, wash gently with water.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C

IV. INFORMATION ABOUT THE WASHING PROCEDURE

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well..
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.
- .

V. PREPARATION OF SAMPLES:

Dilute the samples at **1/50** in a dummy plate (i.e. 5 µl of serum in 245 µl of diluent is enough to run a sample in one well with viral antigen and in another well with control antigen).

VI. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrated washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml of concentrated solution and

960 ml of water). Once prepared this solution remains stable when stored between at +4°C.

- **Control sera (+) and (-):**

Controls are ready to use. Do not dilute. (plate).

- **Conjugate:**

Dilute 1/100 in the conjugate diluent provided in the kit (110µl of concentrated conjugate in 11 ml of diluent is enough for one

VII. TEST PROCEDURE

1. Prior to starting the test, bring all reagents to room temperature (22-25°C).
2. Add 100 µl of the serum samples diluted as previously described onto the positive and negative antigens. **Seal tightly the plate and incubate for 30 min at 37°C.**
3. Wash 3 times following the procedure previously described.
4. Add 100 µL of the conjugate at the indicated dilution, to each well.
5. Wash 6 times following the described procedure.
6. Add 100 µl of the substrate solution. Keep the plate for **10 min at room temperature**. Start timing after adding it to the first well.
7. Add 100 µl of the stop solution, to each well in the same order in which the substrate was added.
8. Read the absorbance of each well with a spectrophotometer at 450 nm within the following 5 min after the addition of the stop solution.

Seal the plate and **incubate for 30 min at 37°C.**

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

The OD obtained with the negative antigen must be subtracted from the OD obtained with the positive antigen (sample and controls)

A. Validation criteria:

The test can be considered valid when:

Positive Control: OD POS Ag – OD NEG Ag > 0,3

Negative Control: OD POS Ag – OD NEG Ag < 0,2

B. Results Interpretation:

- **Positive and negative samples**

Samples should be considered positive for antibodies to BLV when:

OD POS Ag – OD NEG Ag > 0,3

Samples should be considered negative for antibodies to BLV when:

OD POS Ag – OD NEG Ag < 0,25

- **Doubtful samples**

Samples showing values between 0.30 and 0.25 ("Grey zone") should be considered doubtful and be retested. If the result remains the same we recommend to test a new extraction obtained 15 days after the first one.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.com

Distributed in

by: