## INGEZIM BRSV COMPAC

Prod Ref: 12.BRS.K3

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## I. TECHNICAL BASIS

This kit is based on a competitive enzymatic immunoassay. We make a brief description of the technique below:

We fix the BRSV antigen on a solid support (polystyrene plate). Serum samples and a peroxidase conjugated BRSV-specific monoclonal antibody (Mab) are added into wells. If sample contains antibodies against the virus, they will bind to the antigen adsorbed on plate blocking the binding of the conjugated Mab. On the other hand, if the sample does not contain antibodies to the BRSV the conjugated Mab will bind to the antigen.

After washing to eliminate all non-fixed material, we can detect the presence of the conjugated Mab after addition of a substrate-TMB solution. A colorimetric reaction will develop in the case of absence of antibodies in the serum sample and no colour (indicating no binding of the Mab) in case of presence of antibodies against the virus in the assayed sample.

The aim of this kit is to provide the users with a reliable and automatizable diagnostic technique for this disease.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
11. Substrate must be handled with care, it is very sensible to light and contamination.

## III. STORAGE OF COMPONENTS

All reagents and components must be stored between +2°C and +8°C.

**Once opened**, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

We recommend the **TITRATION** of the samples starting from dilution 1/2. Dilutions can be done in the same or in different (blank) plate, keeping in mind that the final volume of the diluted sample in the well

should be 50 µl because you must add 50 µl of conjugate to the same well.

For **SCREENING**, test samples at dilution 1/2.

## VI. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrated washing solution with 24 parts of distilled or deionized water (e. g. 40 ml of concentrate plus 960 ml of water) When ready this solution remains stable when is stored between +2°C and 8°C.

- **Preparation of the control sera:**

We recommend distributing the control sera in aliquots and storing it at -20°C till use.

Control sera needs to be assayed diluted 1/2 with diluent. This dilution could be done directly on each well of the plate by adding 25 µl of serum and 25 µl of diluent.

- **Preparation of the conjugate: to make immediately before use:**

Dilute the needed quantity of conjugate provided in the kit 1/100 into diluent:

→ The necessary and sufficient quantity for a complete plate is 60 µl of conjugate with 6 ml of diluent.

→ The necessary and sufficient quantity for an 8 well strip is 5 µl of conjugate with 0,5ml of diluent.

Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected..

## VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Samples and control sera addition:

**FOR TITRATION:**

Add 50 µl of diluent to each well and 50 µl of samples to the wells of row A (A1.... A12). Then, make two-fold dilutions of each serum up to row H (final dilution is 1/256). Remove 50 µl from the wells of row H. Negative and Positive controls are tested only at 1/2 dilution in duplicate wells (e.g. wells H1, H2 and H3, H4).

**FOR SCREENING:**

Add 25 µl of diluent and 25 µl of sample. Do the same with controls.

3. Add 50 µl of conjugate prepared as previously described to each well. Shake carefully the plate to

help the correct homogenisation of the components. Seal the plate and incubate for **30 min at 37°C**.

4. Wash 6 times following the procedure.
5. Add 100 µl of substrate, to each well. Keep the plate in darkness for **10 min at room temperature**.
6. Add 100 µl of stop solution to each well.
7. Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULTS INTERPRETATION

The reading must be done at **450 nm in an ELISA reader**.

### A. Validation of the test:

The test is considered valid if:

→ The OD for the Positive Control (PC) < 0.3

→ The OD for the Negative Control (NC) range between 0.8 and 1.5

### B. Interpretation of the results:

#### A. Cut Off calculation:

Cut Off = NC - [(NC - PC) x 0.4]

**B. Results Interpretation:**

When you are running duplicate samples, OD values will be calculate for each sample as the arithmetic mean of both values.

→ All samples with OD higher than CUT OFF value must be considered as **negative** to BRSV antibodies.

- All samples with OD values lower than CUT OFF value must be considered as **positive** to BRSV antibodies.
- Samples with OD values ranging  $\pm 5\%$  of cut off value should be considered doubtful and retested.
- Titre of a sample will be the last dilution showing OD value lower than cut off.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: ingenasa@ingenasa.com  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in by: