# INGEZIM BRUCELLA SMALL RUMINANTS

Prod Ref: 13.BM.K1

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## I. TECHNICAL BASIS

Brucella species have two main antigenic components. One of those is a protein-polysaccharide complex LPS-S with endotoxic activity, responsible of the specificity of agglutination and is the major antigenic component. The second one is a polysaccharide (poli-B) with no endotoxic activity and with no relevance role in the agglutination test.

The INGEZIM BRUCELLA kit is based on the indirect immunoenzymatic assay described below:

The antigen is fixed in a solid support (polystyrene plate). The used antigen is a purified extract of the LPS of Brucella Melitensis. The serum sample is added to each well. After incubation and

washing we add a labelled monoclonal antibody (peroxidase conjugated) specific to ruminant immunoglobulins. If the serum sample contains antibodies against the Brucella, the conjugate will bind to the plate, whereas if it does not contain specific antibodies no binding of the conjugate will be appreciated.

After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
11. Substrate must be handled with care, it is very sensible to light and contamination.

## III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

**Once opened**, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

The sera samples need to be at 1/100 dilution (i.e. 5µl of serum with 495 µl of diluent).

## VI. PREPARATION OF REAGENTS

- **Washing solution:**  
Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water (40 ml of concentrated solution and 960 ml of water). When ready this solution remains stable stored at 4°C.
- **Preparation of control sera (+) and (-):**  
Controls are used as samples. You must do the 1/100 dilution in diluent. Once opened, controls are stable between +2°C and +8°C for 4 week. Otherwise we recommended storing the control sera at -20°C in aliquots.
- **Preparation of Diluent:**  
Dilute 1 volume of diluent 5x concentrated supplied in the Kit with 4 volumes of distilled or deionized water (100 ml of concentrate plus 400 ml of water).
- **Preparation of the conjugate: to make immediately before use.**  
Dilute 1/100 with diluent 1x prepared as it is described at point "Preparation of diluent":
  - The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of conjugate in 11 ml of diluent.
  - The necessary and sufficient quantity of conjugate for an 8 wells strip is 10 µl of conjugate in 1 ml of diluent.

Homogenize the solution before use.

**ATTENTION:** Prepare only the quantity needed for each time because the remainder volume has to be rejected.

## VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to bring room temperature before use.
2. Add 100 µl of diluted Positive and Negative control serum to their respective wells,  
  
Add samples (diluted as it's specified in the previous instructions) to remainder wells of the plate. We recommend running samples and control in duplicate for screening purposes. **Incubate the plate for 1 hour at room temperature (20-25°C).**
3. Wash 3 times following the described procedure.
4. Add 100 µl of conjugate (prepared following previous instructions) to each well. **Incubate the plate for 30 minutes at room temperature (+25°C).**
5. Wash 5 times following the described procedure.
6. Add 100 µl of substrate to each well. Keep the plate in darkness for **5 min at room temperature.**
7. Add 100 µl of stop solution to each well.
8. Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

### A. Validation criteria:

The test is considered valid when:

- OD of the positive control is higher than 1.0
- OD of the negative control is lower than 0.25.

### B. Positivity index (IP) calculation:

IP of the sample:

$$\frac{OD_{450} \text{ Sample}}{OD_{450} \text{ C+}} \times 100$$

### C. Results Interpretation:

When you are running duplicate samples, OD values will be calculated for each sample as the arithmetic mean of both values<sup>3</sup>

- All samples with IP ≥ than 25% must be considered as Positives.
- All samples with IP lower than 25% will be considered as Negatives.

Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in \_\_\_\_\_ by: