INGEZIM NEOSPORA 3.0

Prod Ref: 12.NC.K1

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1. Abs C. Pos – Abs C. Neg $\geq 0,4$
2. Abs. del control negativo < 0.25

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I. TECHNICAL BASIS

This enzymatic immunoassay (ELISA) has been developed for the detection of specific antibodies to *Neospora caninum* in bovine serum. Sera samples must be diluted at 1/100 and incubated in the wells coated with *Neospora caninum* antigen for 45 minutes. After the first washing step, the anti-bovine IgG conjugate is added and incubated for 30 minutes. After a second wash, the substrate is dispensed and

after 15 minutes the reaction is stopped. The colour development is read with an ELISA reader. The green colour indicates the presence of specific antibodies to *Neospora caninum*. The more concentration of antibodies it has, the more intensity of colour it develops, and therefore, the optical density (OD) will be higher.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20° - 25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
10. **IMPORTANT!** Substrate must be handled with care, it is very sensitive to light and contamination. Remove the volume to use from the bottle and never put the remaining volume back into it.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic microplate washer or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as it is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing, shake the plate turned over on an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Dilute serum samples 1/100 in diluent DE19-01. This dilution is recommended in two steps as follows: initially dilute 10 µl of serum in 190 µl of diluent and then

transfer 10 µl of this dilution to one well of the ELISA plate contain 40 µl of diluent DE19.

VI. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled water or deionized water (40 ml of concentrate solution and 960 ml of water). This solution remains stable between +2°C and +8°C.

- **Preparation of diluent 1x:**

- Is ready to use. Do not dilute

- **Control sera:**

Control sera are ready to use. Do not dilute

- **Preparation of the conjugate** Is ready to use. Do not dilute

VII. TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use. (Add 50 µl of samples diluted as it is specified in the previous instructions. Add 50 µl of Positive and Negative control serum ready to use to their respective wells. We recommend to run controls in duplicate wells. Cover and incubate the plate for **60 minutes at 37°C ± 1°C**.)
2. Wash 3 times with washing solution following the described procedure.
3. Add 50 µl of Conjugate to each well. Cover and incubate the plate for **60 minutes at 37°C ± 1°C**.
4. Wash 4 times following the described procedure.
5. Add 50 µl of substrate solution, to each well. Shake softly and keep the plate for **20 min at room temperature (23°C ± 2°C)** in a dark place.
6. Add 50 µl of stop solution to each well in the same order as the substrate was dispensed.
7. Read the OD of each well with a spectrophotometer at 450 nm in the following 5 minutes.

VIII. READING AND RESULT INTERPRETATION

B. Positivity Index Calculation (IRPC):

If you run samples in duplicate, OD values will be calculated for each sample as the arithmetic mean of both values.

To obtain IRPC of the sample:

$$\text{IRPC} = \frac{\text{Abs}(\text{Sample}) - \text{Abs}(\text{Control -})}{\text{Abs}(\text{Control +}) - \text{Abs}(\text{Control -})} \times 100$$

C. Validation of the test:

The test could be considered valid if:

1. $\text{OD}_{450} \text{ C (+)} - \text{OD}_{450} \text{ C (-)} \geq 0,4$
2. Negative control $\text{OD} < 0.25$

D. Interpretation of the results:

- Sample $\text{IRPC} \geq 20$ will be considered as positive.
- Sample $\text{S/P} < 17$ will be considered as negative.
- Sample with values between both IRPC will be considered doubtful. A second serum sample could be collected 2 weeks later and retested.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in.....by:

IT-73840
IT-73780

ISO 14001:2015
9191.INGE

ISO 9001:2015
9175.ING2