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
Blocking immunoenzymatic assay for the specific detection and differentiation of antibodies to transmissible Gastroenteritis Virus and Porcine Respiratory Coronavirus in pig serum.

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OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT:

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Distilled or deionised water.
Micropipets from 5 to 200 µl.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 ml
ELISA Reader (450 nm filter)

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I. TECHNICAL BASIS

This kit is based on an enzymatic immunoassay (Blocking ELISA). We make a brief description of the technique below:

The **TGEV recombinant antigen** is captured on the plate through a monoclonal antibody previously fixed on the wells). Each serum sample to assay, needs to be settle in two wells. After incubation we add to one of the wells a specific monoclonal antibody (conjugate A) against one common epitope of both corona viruses. If the serum sample contains antibodies against any one of the viruses, they will not permit the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies the Mab will bind to the antigen on the plate.

In the other well we will add a specific Mab (conjugate B) against one specific

epitope of the TGEV. If the serum sample contains antibodies against this virus (TGEV) the Mab will be competed by these antibodies and it will not bind to the antigen. If serum sample does not contain antibodies against TGEV the Mab will bind to the antigen.

After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled Mabs by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

The different combinations of results in both wells will permit us to know if the serum sample contains antibodies against TGEV, or against PRCoV or do not contain antibodies against any of them.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking or smoking where specimens or Kit reagents are being handle
7. Do not pipette by mouth .Use a new tip for each serum sample.
8. Substrate must be handle with care, it is very sensible to light and contamination.
9. Stop solution is a strong acid. Handle with care.
10. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.

III. STORAGE OF COMPONENTS

Between +2°C and +8°C.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLE

The sera sample needs to be at 1/2 dilution in diluent to be tested.

¡ATTENTION! Take in consideration that sera samples need to be tested in two wells.

For your facility the dilution of sera samples could be done directly on the wells by adding 50 µl of serum and 50 µl of sera diluent to each well.

VI. PREPARATION OF REAGENTS

❖ Washing solution:

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water (i.e. 40 ml of concentrate plus 960 ml of dH₂O). When ready this solution remains stable when is stored between +2°C and +8°C.

❖ Preparation of the control sera:

If you are going to run the kit in a partial way, we recommend distributing the control sera in aliquots and storing it at -20°C till use.

Control sera needs to be assayed at a ½ dilution with serum diluent. This dilution could be done directly on each well of the plate in which the sera will be assayed by adding 50 µl of serum and 50 µl of serum diluent.

**Preparation of the conjugates A and B:
Conjugates are supplied ready to use.
DO NOT DILUTE.**

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in by: