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
Blocking immunoenzymatic assay for antibody detection to Porcine Coronaviruses. This assay does not allow the differentiation of antibodies to PRCV from the TGEV ones.

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OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT:

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Distilled or deionised water.
Micropipets from 5 to 200 µl.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 ml
ELISA Reader (450 nm filter)

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I. TECHNICAL BASIS

This kit is based on a blocking enzymatic immunoassay (Blocking Elisa). We make a brief description of the technique below:

The antigen, a recombinant viral protein is fixed in a solid support (polystyrene plate) by a specific monoclonal antibody. After incubation with serum sample we add a monoclonal antibody (peroxidase conjugated) against one common epitope of both TGE and PRCV viruses. If the serum sample contains antibodies against any one of the viruses, they will not permit the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies the Mab will bind to the antigen on the plate.

After washing the plate to eliminate all non fixed material, we can detect the presence or absence of labelled Mab by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

When the well gives the colorimetric reaction that means that the conjugated Mab has bound to the antigen and for this reason the sample could be marked as negative to Porcine Coronavirus. In the other hand when the well does not give this colorimetric reaction it means that the sample antibodies has bound to the antigen without permitting the binding of the conjugated Mab. In this case the sample must be marked as Positive to Porcine Coronavirus (TGEV and/or PRCV).

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. Do not eat, drink, or smoke where specimens or Kit reagents are being handled.
7. Do not pipette by mouth. Use a new tip for each serum sample.
8. Substrate must be handled with care, it is very sensitive to light and contamination.
9. Stop solution is a strong acid. Handle with care.
10. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.

III. STORAGE OF COMPONENTS

All reagents and components must be stored between +2°C and +8°C.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLE

The sera samples need to be at 1/2 dilution in diluent to be tested.
For your facility the dilution of sera samples could be done directly on the wells by adding 50 µl of diluent and 50 µl of serum to each well.

VI. PREPARATION OF REAGENTS

- Washing solution:**
Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water (i.e. 40 ml of concentrate plus 960 ml of dH₂O).. When ready this solution remains stable when is stored between +2°C and +8°C.
 - Preparation of the control sera (+) and (-):**
If the kit is going to be used partially, we recommend distributing control sera in aliquots and storing it a -20°C till use.
- Control sera must be assayed at a 1/2 dilution with diluent. This dilution could be done directly on each well of the plate in which the sera will be assayed by adding 50 µl of reconstituted serum and 50 µl of diluent.
- Preparation of the conjugates A**
Conjugates are supplied ready to use.
DO NOT DILUTE

VII. TEST PROCEDURE

- All reagents must be allowed to come to room temperature before use.
- Add 50 µl of diluent to each well and 50 µl of Positive control serum, Negative control serum and samples to their respective wells. We recommend running samples and control in duplicate. Incubate the plate **1-hour at 37°C.**
- Wash 3 times.
- Add 100 µl of conjugate to each well.
Incubate the plate for **30 min at 37°C.**
- Wash 4 times following the described procedure.
- Add 100 µl of substrate, to each well. Keep the plate for **10 min at room temperature**, in a dark place.
- Add 100 µl of stop solution to each well.
- Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

A. Validation of the test:

The test could be considered valid if:

The ODs for the Positive Control < 0.3

The ODs for the Negative Control > 1.00

B. Interpretation of the results:

- Cut Off calculation:**

Positive/Negative
Cut Off value = 0.6 x OD of Neg. Control

- Results Interpretation:**

When you are running duplicate samples, OD values will be calculate for each sample as the arithmetic mean of both values.

⇒ All samples with OD higher than CUT OFF value must be considered negative.

- All samples with OD values lower than CUT OFF values must be considered positive

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: