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I. TECHNICAL BASIS

The Kit has been designed to detect, in a very easy way, specific antibodies to TGEV in porcine serum samples. This Kit is based on a blocking enzymatic immunoassay. We make a brief description of the technique below:

The **TGEV recombinant antigen** is captured on the plate through a monoclonal antibody previously fixed on the wells. The serum sample is added to each well. After incubation, we add a specific **monoclonal antibody (peroxidase conjugated) specific to S protein of the TGEV**. If the serum sample does not contain antibodies against the

virus, the conjugate will bind to the plate, whereas if it contains specific antibodies to TGEV, these antibodies will bind the antigen on the plate blocking the binding of the conjugate.

After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

Due to the high specificity of the monoclonal used as conjugate, the assay has a high specificity and sensitivity index.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. Do not eat, drink or smoke where specimens or Kit reagents are being handle
7. Do not pipette by mouth .Use a new tip for each serum sample.
8. Substrate must be handle with care, it is very sensible to light and contamination.
9. Stop solution is a strong acid. Handle with care.
10. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.

III. STORAGE OF COMPONENTS

All reagents and components must be stored between +2°C y +8°C.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLE

The sera sample needs to be at 1/2 dilution in diluent to be tested.

For your facility the dilution of sera samples could be done directly on the wells by adding 50 µl of serum and 50 µl of sera diluent to each well.

VI. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the kit into 24 parts of distilled or deionized water (i.e. 40 ml of concentrate plus 960 ml of dH₂O). When ready this solution remains stable when is stored between +2°C and +8°C.

- **Preparation of the control sera:**

Once opened, control sera are stable for one month. If you are going to run the kit in a partial way, we recommend distributing the

control sera in aliquots and storing it at -20°C till use.

Control sera needs to be assayed at a 1/2 dilution with diluent. This dilution could be done directly on each well of the plate in which the sera will be assayed by adding 50 µl of serum and 50 µl of diluent

- **Preparation of the conjugate:**

Conjugate is provided ready to use.
DO NOT DILUTE

VII. TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. Add 50 µl of diluent to each well.
3. Add 50 µl of PRCV control serum, TGEV control serum and samples to each well of the plate. We recommend running samples and control in duplicate. **Incubate the plate for one hour at 37°C.**
4. Wash 3 times following described procedure
5. Add 100 µl of conjugate to each well. Seal the plate **and incubate for 30 minutes at 37°C**
4. Wash 6 times following described procedure.
5. Add 100 µl of substrate, prepared following previous instruction, to each well. Keep the plate for **10 min at room temperature.**
6. Add 100 µl of stop solution to each well.
7. Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

A) Validation of the test:

The test could be considered valid if:

- The ODs for the TGEV Control Serum < 0.3
- The ODs for the PRCV Control Serum > 0.7

B) Interpretation of the results:

A. Cut Off calculation:

Cut Off (+) = 0.6 × OD of PRCV Control

Cut Off (-) = 0.7 × OD of PRCV Control

B. Results Interpretation:

When duplicate samples are being run, OD values of each sample must be calculated as the arithmetic mean of both values.

Every sample with OD higher than CUT OFF (-) value must be considered as **Negative** to TGEV antibodies.

Every sample with OD value lower that CUT OFF (+) value will be considered as **Positive** to TGEV antibodies.

Samples with OD between these two values should be considered doubtful and retested.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: