INGEZIM BRUCELLA PORCINA

Prod Ref: 11.BP.K.1

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I. TECHNICAL BASIS

Brucella species have two main antigenic components. One of those is a protein-polysaccharide complex LPS-S with endotoxic activity, responsible of the specificity of agglutination and is the major antigenic component. The second one is a polysaccharide (poli-B) with no endotoxic activity and no relevance role in the agglutination test.

The INGEZIM BRUCELA kit is based on the indirect immunoenzymatic assay described below:

The antigen is fixed in a solid support (polystyrene plate). The used antigen is a

purified extract of the LPS of Brucella. The serum sample is added to each well. After incubation and washing we add a labelled monoclonal antibody (peroxidase conjugated) specific to pig IgG. If the serum sample contains antibodies against the Brucella, the conjugate will bind to the plate, whereas if it does not contain specific antibodies no binding of the conjugate will be appreciated.

After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the substrate that in presence of the peroxidase will develop a colorimetric reaction.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Substrate and substrate buffer must be handled with care.
10. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
11. Stop solution may precipitate at +4°C. When it happens, warm at 37°C before use.

III. STORAGE OF COMPONENTS

Plates and reagents must be stored between +2°C and +8°C.

Once opened, control sera are stable for one month. In case that they are not going to be used in this period, we recommend to store them at -20°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

INGEZIM BRUCELLA PORCINA 11.BP.K.1

V. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water. Once ready, this solution remains stable at +4°C.

- **Preparation of the conjugate (to make immediately before use):**

Dilute the needed quantity of conjugate provided in the Kit 1/100 with diluent:

- The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of conjugate in 11 ml of diluent.

- The necessary and sufficient quantity of Conjugate for an eight wells strip is 10 µl of conjugate in 1 ml of diluent.

Mix thoroughly the solution before use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

- **Control Sera:**

Controls are used as samples. You must do the 1/25 dilution in diluent. Once opened, controls are stable between +2°C and +8°C for 4 week. Otherwise we recommended storing the control sera at -20°C in aliquots.

VI. PREPARATION OF SAMPLES

The sera samples need to be at 1/25 dilution. We recommend to do this dilution in an empty plate (i.e. 240µl of diluent and 10 µl of each serum, samples or controls). Shake well in order to homogenate and add 100 µl / well.

VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 100 µl of diluted Positive and Negative control serum to their respective wells. Add samples (diluted as it's specified in the previous instructions) to remainder wells of the plate. We recommend (not mandatory) running samples and control in duplicate for confirmatory purposes. **Incubate the plate for 1 hour at room temperature (20-25°C).**
3. Wash 3 times following the described procedure.
4. Add 100 µl of conjugate (prepared following previous instructions) to each well. **Incubate the plate for 30 minutes at room temperature (+25°C).**
5. Wash 5 times following the described procedure.
5. Add 100 µl of substrate to each well. Keep the plate in darkness for **5 min at room temperature.**
6. Add 100 µl of stop solution to each well.
7. Read the OD of each well with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

A. Validation of the test:

The test is considered valid when:

OD value of the positive control is higher than 1.0 and the OD value of the negative control is lower than 0.3

B. S/P calculation:

$$S/P = DO \text{ sample} / DO \text{ control} +$$

C. Interpretation of the results:

When you are running duplicate samples, OD values will be calculated for each sample as the arithmetic mean of both values

- All samples with S/P higher than or equal to 0.25 should be considered as Positive.
- All samples with S/P value lower than 0.25 should be considered as Negative.
- LPS of Brucella and Yersinia O9 share similar antigenic determinants so, some cross reaction could be obtained.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91.368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: