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
$$\text{TITULO (Y)} = 156 (e^{3,6 \times 0,333}) = 517$$

I. TECHNICAL BASIS

The Kit has been designed to detect, in a very easy way, specific antibodies to E. rhusiopathiae in pig

sera samples. This Kit is based on an indirect enzymatic immunoassay. We make a brief description of the technique bellow:

The antigen is fixed in a solid support (polystyrene plate). The serum sample is added to each well. After incubation we add a monoclonal antibody (peroxidase conjugated) specific to porcine immunoglobulins. If the serum sample contains antibodies against the bacteria, the conjugate will bind to the plate, whereas if it does not contain specific antibodies no binding of the conjugate will be appreciated.

After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the substrate that in

presence of the peroxidase will develop a colorimetric reaction.

The major advantage of this assay is that could be applied not only for screening proposes, but also for titration. By this way seroconversion process of the animals could be monitored easily.

The aim of this kit is provide to the users a reliable and automatizable diagnostic technique for this disease that till this time is been made by the agglutination test.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit use a fresh substrate preparation.
10. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
11. Substrate and substrate buffer must by handle whit care.

III. STORAGE OF COMPONENTS

Store all plates and reagents between +2°C y +8°C.

ATTENTION! Stop solution may precipitate at +this temperature, if it happens, warm at 37°C till total dissolution.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichanel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

For **screening and/or one well titration**, serum samples must be tested at dilution 1/200 in diluent.

E.I. 2 µl of serum sample in 400 µl of diluent.

For **titration purposes**, we recommend to assay 4 two-fold dilutions for each sample (1/200,

1/400, 1/800 and 1/1600). For your facility, this step could be done directly on the plate.

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water. When ready this solution remains stable between +2°C and +8°C.
 - **Diluent:**
Dilute 1 volume of diluent 5x concentrated supplied in the Kit with 4 volume of distilled or deionized water (i.e. 100 ml of diluent with 400 ml of water). Store this solution between +2°C and +8°C.
 - **Preparation of the conjugate: to make immediately before use:**
Dilute the needed quantity of conjugate provided in the Kit 1/100 with diluent:
 - ⇒ The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of Conjugate with 11 ml of diluent.
 - ⇒ The necessary and sufficient quantity of Conjugate for an eight wells strip is 10 µl of conjugate with 1 ml of diluent.
- Shake very well the solution before use.
Prepare only the quantity needed for each time because the remainder volume has to be rejected.
- **Control Sera:**
Controls are supplied ready to use. For titration make two fold dilution
 - **Substrate:**
Ready to use

VII. TEST PROCEDURE:

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 100 µl of Positive and Negative control sera, to two wells. Add 100 µl of each dilution of serum samples to be assayed to the remainder wells. **Seal the plate and incubate 1 hour at 37°C.**
3. Wash 3 times following the described procedure.
4. Add 100 µl of conjugate (prepared following previous instructions) to each well. Shake. **Seal the plate and incubate for 45 at +37°C.**
5. Wash 3 times following the described procedure.
6. Add 100 µl of substrate solution to each well. Keep the plate for **15 min at room temperature.**
7. Add 100 µl of stop solution to each well
8. Read at 405 nm within 5 min after the addition of stop solution.

VIII. READING AND INTERPRETATION OF THE RESULTS:

The reading must be done with a spectrophotometer at 405 nm.

a). Validation of the assay:

The assay will be valid only when:

1. OD of positive control is higher than 1,500.
2. OD of negative control is lower than 0,300.

b). Cut off calculation:

- ⇒ Positive Cut off = OD of negative control plus 0,200
- ⇒ Negative Cut off = OD of negative control plus 0,100.

c). Result interpretation:

- For screening purposes:
 - All samples with an OD value higher than positive cut off must be considered positives.
 - All samples with an OD value lower than negative cut off, must be considered negatives.
 - Samples with OD values between both Cut off values, must be considered as doubtful. In this cases a new test of the samples is needed.
- For titration purposes:
The titre of each sample will be the highest dilution which develop a positive result (OD value higher that Cut off value).

d). Example

- 1). Assay date:
 Mean OD Control + = 2.130
 Mean OD Control - = 0.105
 OD Sample 1 1/200 = 0.709
 OD Sample 1 1/400 = 0.372
 OD Sample 1 1/800 = 0.242
 OD Sample 1 1/1600 = 0.164
- 2 Validation of the assay
 OD Control + = 2.130 (> 1.5) OK
 OD Control - = 0.102 (< 0.3) OK
- 3). Cut off calculation
 Positive cut off = 0.105 + 0.2 = 0.305
 Negative cut off = 0.105 + 0.1 = 0.205
- 4). Interpretation
 Sample 1: Positive (0.709)
 Titer Sample 1 : 1/400 (0.372)

IX. ONE WELL TITRATION:

Titre of positive samples will be obtained applying the following formula:

$$\text{TITRE (Y)} = 156 (e^{3,6x})$$

Where **e** is the natural logarithm base (2.718282); **x** is the sample's OD expressed as the relationship S/P (sample's OD / Positive control OD)

EXAMPLE:

3. Results of the assay:

Mean of Positive control OD = 2.130
 Mean of Negative control OD = 0.105
 Mean of OD for Sample 1 at dilution 1/200 = 0.709
 S/P = 0,709/2,130 = 0,333

4. Interpretation:

Sample 1 is positive because the OD at dilution 1/200 is higher than the positive cut off. (see example of point VIII)
 The titre of the sample is:

$$\text{TITRE (Y)} = 156 (e^{3,6*0,333}) = 517$$

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: