INGEZIM ARTERITIS 2.0

Prod Ref: 14.EA2.K1

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I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). The technique is briefly described below:

We fix positive and negative antigens on different plates. When a sample serum contains specific antibodies against the virus, they will bind the viral antigen adsorbed on plate. After washing the plate to eliminate all non-fixed material from the sera sample, the presence of horse immunoglobulins can be detected

using a specific peroxidase conjugated. After the addition of the substrate a colorimetric reaction will appear that can be measured using an ELISA reader.

Therefore the difference between the OD obtained with the positive and the negative antigen will determine the presence or absence of specific antibodies against the virus in the horse serum.

The aim of this kit is to provide users with a reliable and fast diagnostic technique for this disease, which up until now has been made by virus neutralization.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix reagents or use instructions from different kits.
4. Avoid any contamination of the kit reagents.
5. Do not use components after expiry dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking while handling reagents or samples.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. The stop solution is a strong acid and therefore must also be handled with care. In case of an accidental contact with skin, rinse gently with water.
10. Systematically include a positive control and a negative serum each time the assay is run.
11. The substrate must be handled with care, as it is very sensitive to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps can be done using an automatic plate washer or a multichannel pipetting device suitable for dispensing 300 µl in each well.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over, in order to avoid the possible exchange of contents from one well to another.
- Dispense a volume of 300 µl of washing solution in each well.
- Shake the plate gently, avoiding cross-contamination between wells.
- Briskly turn the plate over to empty the wells.
- Repeat the process as many times as indicated in the kit instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to be used. Do not let the plate dry longer than strictly necessary.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

V. PREPARATION OF SAMPLES:

Sera samples must be assayed at 1/100 dilution. (i.e. 5 µl of serum in 495 µl of diluent is enough to run a sample in duplicate)

VI. PREPARATION OF REAGENTS

- **Washing solution:**

Dilute one part of the concentrated washing solution provided in the Kit with 24 parts of distilled or deionized water (40 ml of concentrate solution and 960 ml of water). Once prepared, this solution remains stable when stored between +2°C y +8°C.

- **Control sera (+) and (-):**

Controls are ready to use. Do not dilute.

- **Conjugate::**

Controls are ready to use. Do not dilute.

VII. TEST PROCEDURE

1. Prior to starting the test, bring all reagents to room temperature (20-25°C) (except conjugate).
2. Add 100 µl of the sera samples, diluted as previously described, to each of the appropriate wells in both the positive and negative plate. Then add 100 µl of the positive and negative control sera to both plates in their corresponding wells. **Seal tightly the plate** (to avoid evaporation) and **incubate for 2 h at 37°C**.
3. Wash 3 times following the procedure previously described.
4. Add 100 µl of the conjugate to each well.

Seal the plate and **incubate for 30 min at room temperature**.
5. Wash 6 times following the procedure previously described.
6. Add 100 µl of the substrate solution. Keep the plate for **10 min at room temperature**. Start timing after filling the first well.
7. Add 100 µl of the stop solution, to each well following the same order in which the substrate was added.
8. Read the OD of each well at 450 nm within the following 5 min after the addition of the stop solution 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

The OD obtained with the negative antigen must be subtracted from the OD obtained with the positive antigen (for samples and controls)

A. Validation criteria:

The test can be considered valid when :

Positive Control: OD POS Ag – OD NEG Ag > 0,5

Negative Control: OD POS Ag – OD NEG Ag < 0,1

B. Results Interpretation:

- **Positive and negative samples**

A sample should be considered positive for antibodies to EAV when

OD POS Ag – OD NEG Ag > 0,15

A sample should be considered negative for antibodies to EAV when

OD POS Ag – OD NEG Ag < 0,1

- **Doubtful samples**

Samples showing values between 0.1 and 0.15 ("Grey zone") should be considered doubtful and be retested. If the same result is obtained we recommend assaying a new extraction taken 15 days after the first one.

All positive and doubtful (grey zone) samples must be confirmed by seroneutralization.

Desarrollado y fabricado en España por:
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: