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INGEZIM MIXOMATOSIS 17.MIX.K1

I. TECHNICAL BASIS

The kit is based on the indirect immunoenzymatic assay technique described below:

The antigen is fixed in a solid support (polystyrene plate). The serum sample is added to each well. After incubation and a washing step, a peroxidase conjugate is added. If serum sample contains antibodies specific of myxomatosis virus, the conjugate will bind to them,

whereas if it does not contain specific antibodies no binding will be appreciated.

After washing the plate to eliminate all non fixed material, we can detect presence or absence of conjugate by adding the substrate which, in presence of the peroxidase, will develop a colorimetric reaction. The kit could be used both for screening and for titration purposes.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
10. **IMPORTANT** TMB Substrate must be handle with care, it is very sensible to light and contamination: pipette a sufficient amount for the assay from the substrate storage bottle into a separate dark container prior to colour develop step.
11. Stop solution is a strong acid. Handle with care

III. STORAGE OF COMPONENTS

Store all plates and reagents between +2 and +8°C,

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichanel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Before assaying, sera samples must be diluted 1/200 in the supplied Diluent (5 µl of serum sample to 1 ml with diluent 1X).

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water(40 ml of concentrate solution and 960 ml of water). Once this solution ready, it remains stable between +2°C and +8°C.
- **Diluent:**
Dilute one part of the concentrate diluent provided

in the Kit with 4 parts of distilled or deionized water.

- **Control Sera**
Controls are ready to use.
- **Conjugate:** Ready to use. Do not dilute

VII. TEST PROCEDURE

1. Take out the plate (or strips to be used) from the refrigerator and bring to room temperature before starting the test.
2. Add 100 µl of positive control serum and 100 µl of negative control serum in duplicate wells. Add 100 µl of each sera sample (prepared as specified) to each remainder wells of the assay plate. Seal the plate and **incubate for 1 hour at 37°C**.
3. Wash 3 times following previous instructions
4. Add 100 µl of conjugate to each well. Seal the plate and **incubate 1 hour at 25°C (room temperature)**.
5. Wash 5 times following the described procedure.
6. Add 100 µl of substrate, to each well. Keep the plate **for 10 min at room temperature**.
7. Add 100 µl of stop solution to each well. It is recommended to add this reagent in the same order than was used for adding the substrate.
8. Read the OD of each well at 450 nm within 5 min. after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done at 450nm. If you run controls and samples in duplicate wells, the arithmetic mean of the OD must be considered.

A. Validation of the Test:

OD Positive Control / OD Negative Control > 5.

B. Interpretation of the Results:

A. Cut off calculation:

Cut off Negative

OD Negative control + 0.15

Cut off Positive

OD Negative control + 0.2

All samples with OD higher than cut off positive must be considered as positives, and samples with OD lower than cut off negative, must be considered as negatives. Samples with OD between both values must be considered doubtful.

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If you want to know the titre of positive samples:

- 1) Calculate the ratio OD Sample/OD control + (S/P)
- 2) Substitute "x" in the formula with this S/P value
Titre = 2000 . x^{1.4}

EXAMPLE

Positive Control OD = 1.69

Negative Control OD = 0.15

Cut off Negative = 0.15 + 0.15 = 0.300

Cut off Positive = 0.15 + 0.2 = 0.350

Abs **Sample A** = 1.02 (>0,35 , positive)

Titre = 2000 x 0,603^{1.4} =985

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: