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
I. TECHNICAL BASIS

The kit have been designed to detect antibodies specific for LPS of Brucella spp. being able to detect a very low titters of antibodies in sera of different infected animal species (bovine, ovine, porcine, etc).

The INGEZIM BB COMPAC kit is based on the blocking enzyme immunoassay described bellow:

The solid phase is plates coated with purified LPS of Brucella abortus. After adding the sample to the well, if it contains specific antibodies against Brucella, they will bind to the antigen absorbed on plate while if the sample does not contain specific antibodies they will not bind to the antigen.

If we add a specific monoclonal antibody (conjugated with peroxidase) against LPS antigen coated to the plate, it will compete with the antibodies of the serum. If the serum samples contain specific antibodies, they will not permit the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies the Mab will bind to the antigen on the plate. After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled Mab by adding the substrate (TMB) that in presence of the peroxidase will develop a colorimetric reaction.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative serum must be tested in a systematic way.
11. Substrate must be handle whit care, it is very sensible to light and contamination.

III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichanel pipetting device suitable for dispensing 300 µl on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

- a) Small ruminants (ovine, caprine):** Sera samples must be tested at 1/5 dilution in diluent. This dilution could be made directly on the assay plate by adding 80 µl of diluent and 20 µl of sample to each well.
- b) Bovine and Porcine:** Samples must be tested at 1/10 dilution (90 µl of diluent and 10 µl of serum). Mix gently to homogenate the dilution in the well.

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionised water (40 ml of the concentrated solution in 960 ml of distilled or deionised water). Once ready this solution remains stable between +2°C and +8°C.
- **Control sera and conjugate**
They are supplied ready to use. DO NOT DILUTE.

VII. TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. Add samples as indicated in point V. Add 100 µl of controls in duplicate wells. Gently mix the content of the wells, seal the plate and incubate for **1 h at room temperature (20-25°C)**
3. Wash 4 times following the described procedure.
4. Add 100 µl of conjugate to each well. Seal the plate and incubate for **1 hour at room temperature (20-25°C)**
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate, to each well. Keep the plate for **10 min at room temperature.**
7. Add 100 µl of stop solution to each well.
8. Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

In case that the samples had been run in duplicate, it has to be considered the mean of both OD values. In the same way, you have to make the mean of the values obtained in the two wells of positive control and the two wells of negative control.

A. VALIDATION OF THE RESULTS

- The OD in the wells with Negative Control (NC) must be higher than 1.
- The OD in the wells with Positive Control (PC) must be lower than 0.35.

B. INTERPRETATION OF THE RESULTS

Calculate the Percentage of inhibition (PI) of each sample as follow:

$$PI = 100 \times [1 - (OD \text{ sample} / OD \text{ negative control})]$$

Samples with PI ≥ 40% should be considered positive for brucella antibodies

Samples with PI < 40% should be considered negative for Brucella antibodies

IX. ADDITIONAL INFORMATION:

This assay can be used as it is described below:

- ✓ Porcine and Bovine species: technique approved for official diagnostic. specie bovina y porcina: Técnica aprobada para el diagnóstico oficial según normativa (Directive 64/432/CEE, Decision 2012/137/UE).
- ✓ Ovine and caprine species: Use for non vaccinated animals in free areas. This technique has not been approved yet for official diagnostic by OIE. Only for research use.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tif: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: