INGEZIM CORONA FELINO

Prod Ref: 16.FCV.K1

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I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the antigen on a polystyrene plate. When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on the plate. After washing to eliminate all non-fixed material from the sera sample, we can detect the presence of cat antibodies using a specific peroxidase conjugate.

After addition of the substrate a colorimetric reaction will appear which could be measured by a spectrophotometer.

In this way the presence of colour means the presence of antibodies against the virus in the cat sera, and the absence of colour the absence of specific antibodies.

The aim of this kit is provide to the users with a reliable and plain diagnostic technique for this disease.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **IMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.
9. Stop solution is a strong acid. Handle with care.
10. Substrate must be handled with care, it is very sensible to light and contamination.

III. STORAGE

Storage: All reagents and plates must be stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Samples must be cat sera. Make a 1/200 dilution (5 µl of serum in 1 ml of diluent). Add 200 µl of this dilution on the first well of a row and then make two-fold dilutions by passing 100 µl to the next wells (previously filled with 100 µl of diluent).

VI. PREPARATION OF REAGENTS

- **Washing solution:**
Dilute one part of the concentrate washing solution provided in the kit with 9 parts of distilled or deionized water. When ready this solution remains stable at +4°C.
- **Control sera:**
Controls are ready to use.
- **Preparation of the conjugate: to make immediately before use.**
Dilute the needed quantity of conjugate provided in the kit 1/100 into the supplied diluent (e.g. 10 µl of conjugate into 990 µl of diluent is enough quantity for one strip).

Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected..

VII. TEST PROCEDURE:

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add samples diluted as previously described to each well to be used. Add positive and negative controls to the wells reserved for them. Seal the plate and **incubate for 30 minutes at room temperature** (20 - 25°C).
3. Wash 4 times following the described procedure.
4. Add 100 µl of conjugate, prepared following the 8. previous instructions, to each well. Seal the plate and incubate for 15 **minutes at room temperature** (20-25°C).
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate solution to each well. Keep the plate for **5 minutes at room temperature**.
7. Add 100 µl of stop solution to each well.

Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

Kit Validation:

The test is validated when:

OD (Negative Control): ≤ 0.20
OD (Positive Control)= 0.9 ± 0.3

Results Interpretation:

Taking the negative control OD value, will be determined the following cut off:

Cut off:

OD of Negative Control x 4

Coronavirus-Positive samples:

Samples with an OD higher than the Cut Off value.

Coronavirus-Negative samples:

Samples with an OD lower than the Cut Off value.

Titre of the samples will be the last dilution showing an OD higher than the cut off.

75% of the FIPV-infected cats show titres $\geq 1/6400$.
50% of the FIPV-infected cats show titres $\geq 1/12800$.

On the other hand, only 20% of FECV-infected cats show titres of 1/6400 and only 5% of them, titres of 1/12800.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA
APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

	Distributed in by:
--	-----------------------