## INgezim FIV IC

Prod. Ref: **1.6.FI2.K41**

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## INgezim FIV IC

Prod. Ref: **1.6.FI2.K41**

Immunocromatographic device for detection of specific antibodies to feline Immunodeficiency virus in feline blood, serum or plasma samples.

INMUNOLOGIA Y GENETICA APLICADA, SA (INGENASA)

C/ Hnos Garcia Noblejas, 39 – 28037 MADRID (Spain)

Tel: (+34)91 3680501

www.ingenasa.com

9001: IT-73780 ISO 14001 ISO 9001  
14001: IT-73840 9191.INGE 9175.ING2

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### KIT COMPONENTS:

| Componente                       | Unidades |    |
|----------------------------------|----------|----|
|                                  |          |    |
| Immunocromatographic devices     | 30       | 12 |
| Plastic disposable pipettes      | 30       | 12 |
| Droppers containing test diluent | 2        | 1  |

### OTHER MATERIALS OR REAGENSTS NEEDED NOT PROVIDED:

– None

### KIT CONSERVATION:

All components of the kit should be stored between + 4 ° C and + 25 ° C, in its original packaging until use.

### PRINCIPLE OF THE TEST:

The test is based on the direct migration technique, and has been designed to detect specific antibodies to FIV in blood, serum and plasma feline samples.

### PRECAUTIONS:

- Carefully follow the instructions.
- Place all reagents to room temperature at least 20 min before use.
- Never open the foil pouch containing the device if not used immediately. It could be hydrated and become unusable. Open only immediately before use.
- Do not mix instructions or reagents from different kits.
- Avoid any contamination of reagents.
- Do not use the kit once exceeded their expiration date.
- Do not eat, drink or smoke when samples and / or reagents are handled.
- Do not touch the chromatographic membrane with any object, especially with the dropper tip.

### INSTRUCTIONS:

#### A.- SAMPLE PREPARATION:

The sample of choice for diagnostic purposes, shall be complete blood, serum or plasma.

They should be fresh and not cloudy. Serum samples may be stored in the refrigerator between 2-8°C for 1-2 days. For longer storage, samples should be stored frozen at -20°C. If working with frozen samples, they should be previously thawed, brought to room temperature and homogenized before their addition to the test.

#### B.- TEST PROCEDURE:

Remove the diagnostic device from its envelope and place it on a flat surface. Discard the bag of silica gel. NOTE: Only open the bag at the time of going to test.

Using the provided pipette, add 1 drop of the sample into the window "S".


**CAUTION: The pipette should be placed vertically on the window and keeping the tip 1 cm to avoid touching the membrane.**

Add 2 drops of test diluent into the window "S".

#### C.- READING AND RESULTS INTERPRETATION:

Keep the device on a flat surface and untouched for 10 minutes, after which proceed to reading the result:

It is advisable not to give any negative sample definitely until you have 10 minutes elapsed from the placement of the sample into the device.