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## I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the antigen on a polystyrene plate. When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on the plate. After washing to eliminate all non-fixed material from the sera sample, we can detect the presence of cat antibodies using a specific peroxidase conjugate.

After addition of the substrate a colorimetric reaction will appear which could be measured

by an spectrophotometer.

In this way the presence of colour means the presence of antibodies against the virus in the cat sera, and the absence of colour the absence of specific antibodies.

The aim of this kit is provide to the users with a reliable and plain diagnostic technique for this disease.

In our kit is important to remark the use of a synthetic peptide like antigen. This method warranties the total absence of infectivity in the kit.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use. **IMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.
9. Stop solution is a strong acid. Handle with care.
10. Substrate must be handled with care, it is very sensible to light and contamination.

## III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

The sera samples need to be at **1/100 dilution** in diluent to be tested.

For titration we recommend assaying different two fold dilutions, starting from 1/100.

## VI. PREPARATION OF REAGENTS

### Σ WASHING SOLUTION:

Dilute one part of the concentrate washing solution provided in the kit into 9 parts of distilled or deionized water (i.e. 100 ml of concentrate and 900ml of water). When ready this solution remains stable when is stored between +2°C and +8°C.

### Σ PREPARATION OF THE CONJUGATE: TO MAKE IMMEDIATELY BEFORE USE.

Dilute 1/100 with diluent.

→ The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of conjugate in 11 ml of diluent.

→ The necessary and sufficient quantity of conjugate for an 8 wells strip is 10µl of conjugate in 1 ml of diluent.

Prepare only the quantity needed for each time because the remainder volume has to be rejected.

Shake very well the solution prior to use.

## VII. TEST PROCEDURE:

- All reagents (except conjugate) must be allowed to come to room temperature before use.
- Add 100 µl of positive and negative controls on two wells. Add 100 µl of each dilution of sera samples to test, prepared following the previous instructions, on the remainder wells of the plate (it is important to run duplicates for both controls and samples test). Seal the plate and incubate for **10 minutes at room temperature.**
- Wash 4 times following the procedure previously described.
- Add 100 µl of conjugate, prepared following the previous instructions, to each well. Seal the plate and incubate for **10 minutes at room temperature.**
- Wash 4 times following the described procedure.
- Add 100 µl of substrate solution, to each well. Keep the plate for **5 min** at room temperature.
- Add 100 µl of stop solution to each well.
- Read the absorbances of each well with an spectrophotometer at 450 nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**.

Samples with an OD higher than the positive cut off value.

### Kit Validation:

The test could be validated when:

- O.D 450 nm positive control >1.**
- O.D 450 nm negative control <0,2.**

### NEGATIVE SAMPLES:

Samples with an OD lower than the negative cut off value.

### Results Interpretation:

Taking the negative control value, the positive and negative cut off values must be calculated as follows:

$$\frac{\text{NEGATIVE CUT OFF:}}{\text{O.D 450 nm negative control} + 0,25}$$

$$\frac{\text{POSITIVE CUT OFF:}}{\text{O.D 450 nm negative control} + 0,3}$$

### DOUBTFUL SAMPLES:

Samples with an OD between both cut off. In these occasions, it is suggested to repeat the test after 3-4 weeks. If there is no seroconversion and its OD is not higher than the positive cut off in this second test, the sample will be considered negative.

If you run the samples in duplicate wells, the OD value for the sample will be the mean of both wells.

### POSITIVE SAMPLES:

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tif: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in \_\_\_\_\_ by: