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
## I. TECHNICAL BASIS

---

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the specific antigen on a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on plate. After washing to eliminate all non fixed material from the sera sample, we can detect the presence of dog immunoglobulins using an specific peroxidase conjugate.

After addition of the substrate a colorimetric reaction will appear which could be measured by an spectrophotometer.

In this way the presence of colour means the presence of antibodies against the virus in the dog sera , and the absence of colour the absence of specific antibodies.

The aim of this kit is provide to the users with a reliable and plain diagnostic technique for this disease.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

---

1. Read the instructions of use carefully.
2. **Bring all reagents to room temperature (20°-25°C) prior to use. iIMPORTANT!**
3. Do not mix instructions or reagents from different kits.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
6. Do not pipette by mouth.
7. Use a new tip for each serum sample.
8. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.
9. Stop solution is a strong acid. Handle with care.
10. Avoid any contamination of the reagents of the Kit. Take special care of substrate because is very sensitive to light and contamination.

## III. STORAGE OF COMPONENTS

---

All reagents and plates must be stored between +2°C y +8°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

---

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

The sera samples need to be at **1/100 dilution** in diluent to be tested.

## VI. PREPARATION OF REAGENTS

- **Washing solution:**  
Dilute one part of the concentrate washing solution provided in the kit into 9 parts of distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

## VII. TEST PROCEDURE:

1. **ALL REAGENTS AND WELLS TO BE USED MUST BE ALLOWED TO COME TO ROOM TEMPERATURE BEFORE USE (AT LEAST FOR 1 HOUR).**
2. Add 100 µl of positive and negative controls to two different wells. Add 100 µl of each dilution of sera samples to test, prepared following the previous instructions, on the remainder wells of the plate. We recommend running the test (sample and control sera), in duplicates. Seal the plate and **incubate for 10 minutes at room temperature.**
3. Wash 4 times following the procedure previously described.
4. Add 100 µl of conjugate to each well. Seal the plate and **incubate for 10 minutes at room temperature.**
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate solution, to each well. Keep the plate for **5 min at room temperature.**
7. Add 100 µl of stop solution to each well.
8. Read the absorbances of each well with an ELISA reader at 450 nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading must be done with an ELISA reader at **450 nm**.

### Kit Validation:

The test is valid when:

- **O.D positive control is >1**
- **O.D negative control is < 0,2**

### Results Interpretation:

Using the positive control value, the cut off points values must be calculated as follows:

#### CUT OFF:

$$\text{O.D 450 nm positive control} \times 0.2$$

If you run the samples in duplicate wells, the OD value for the sample will be the mean of both wells.

#### + **NEGATIVE SAMPLES:**

Samples with an OD lower than the Cut Off value

#### + **POSITIVE SAMPLES:**

Samples with an OD higher than the Cut Off value. In this case you could distinguish three kind of samples

- **low titter** (corresponding to IFI values of 1/20-1/40). These samples show OD values between:  
0.2 x OD of positive control and  
0.4 x OD of positive control.
- **medium titres** (corresponding to IFI values of 1/80-1/160). These samples show OD values between:  
0.4 X OD of positive control and  
0.8 x OD of positive control
- **high titres** (corresponding to IFI values  $\geq$  1/320). These samples show OD values higher than,  
0.8 x OD of positive control

**In spite of the good sensitivity and specificity of this test, all the positive samples should be confirmed with other diagnostic method. Definitive diagnostic must not be based on a unique test result. These results must be correlated with other data: animal history, vaccination and clinical symptoms.**

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: ingenasa@ingenasa.com  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in

by: