INGEZIM PARVOVIRUS IgM

Prod Ref: 15.CPM.K2

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I. TECHNICAL BASIS

This kit is based on a capture enzymatic immunoassay. The technique is briefly described below:

If both, primary vaccinations and infections, a high rate of IgM antibodies to CPV are produced. These are the antibodies which are detected in our assay. Plates are coated with a monoclonal antibody specific to dog IgM.

Sera samples are placed in each well. The existing IgM antibodies present in the serum will be captured by the mAb which is coating the plate. After a washing step to remove any unbound material, a recombinant viral antigen should be added and will remain fixed to the well only if the

sera sample contained specific IgMs to canine Parvovirus.

After additional washing steps in order to remove again any unbound material, the presence of the recombinant antigen can be proved by adding a mAb conjugate specific to PPV labelled with peroxidase.

With the further addition of an adequate substrate, those wells with sera samples with containing specific IgM, will show a colorimetric reaction, on proportion to the titre or level of antibodies, whilst on the other hand, those wells where negative sera samples have been ran, no colorimetric reaction will be observed.

II. PRECAUTIONS AND WARNINGS FOR USERS:

Read the instructions of use carefully.

1. Bring all reagents to room temperature (20°-25°C) prior to use. **IMPORTANT!**
2. Do not mix reagents or use instructions from different kits.
3. Do not use components after expiration dates and do not mix components from different batches.
4. There should be no eating, drinking, or smoking while handling reagents and/or samples.
5. Do not pipette by mouth.
6. Use a new tip for each serum sample.
7. Include a positive and negative control each time the assay is run.
8. The stop solution is a strong acid and should be handled with care.
9. The substrate is very sensitive to light and contamination. It is recommended to pipette a sufficient amount for the assay from the substrate storage bottle into a separate container and discard the remaining volume (never return it to the bottle!).

III. STORAGE OF COMPONENTS

All reagents supplied with the kit, should be kept refrigerated between +2 ° C and +8 ° C until use.

IV. INFORMATION ABOUT THE WASHING PROCEDURE

The washing of the plates can be done using an automatic plate washer, a multichannel pipetting device suitable for dispensing 300 µL on each well, or a squeeze bottle.

After the incubation periods, the washing steps must be done following the instructions below:

- Throw away the content of the plate by briskly turning the plate over, in order to avoid the possible exchange of contents from one well to another.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake the plate gently, avoiding cross-contamination between wells.
- Turn over the plate briskly to empty the wells.
- Repeat the process as many times as indicated in the kit's instructions.
- Prior to emptying the content of the plate in the last washing step, verify that the next reagent to be added is ready to be used. Do not let the plate dry longer than strictly necessary.
- After the last washing step tap the plate upside down on absorbent filter paper to remove any remaining washing solution.

V. PREPARATION OF SAMPLES:

Sera samples must be tested at 1/100 dilution in serum diluent (eg. 5µl of serum in 500µl. of diluent)

INGEZIM PARVOVIRUS IGM 15.CPM.K.2

VI. PREPARATION OF REAGENTS

- **Washing solution :**
Dilute one part of the concentrated washing solution provided in the kit into 9 parts of distilled or deionized water. Once prepared this solution remains stable when stored between +2°C and +8°C.
- **Antigen:**
Before its addition to the well, dilute 1/10 in the diluent provided. (1.1 ml of the recombinant antigen with 10 ml of diluent is enough for one plate. For an 8 well strip dilute 0.1 ml of antigen in 0.9 ml of diluent).
- **Conjugate**
Dilute 1/100 in the conjugate diluent. It is advisable to dilute only the amount of conjugate necessary to run the test, as the remaining solution must be discarded: (110 µl of conjugate in 11 ml of diluent is sufficient for one plate. For an 8 well strip, dilute 10 µl of conjugate in 1 ml of diluent).
- **Controls**
Controls are ready to use: DO NOT DILUTE

VII. TEST PROCEDURE:

1. Prior to starting the test, bring all reagents to room temperature (22-25°C). (except for the conjugate and the recombinant antigen)
2. Add 100 µL of each sample (prepared according to the instructions above). Then add 100 µL of the positive and negative controls. In order to obtain reliable results, we recommend running samples and controls in duplicate wells. Seal the plate and incubate 15 minutes at 37°C.
3. Wash 4 times following the washing procedure, previously described.
4. Add 100 µL of the antigen, prepared according to previous instructions, to each well. Seal the plate and incubate 15 minutes at 37°C).
5. Wash 4 times as indicated before.
6. Add 100 µL of the conjugate at the indicated dilution, to each well. Seal the plate and incubate 15 minutes at 37°C.
7. Wash 4 times.
8. Add 100 µL of the TMB substrate to each well. In order to speed up this process, it is advisable to use a multichannel pipette. Keep the plate for 5 min at room temperature.
9. Add 100 µL of the stop solution to each well. We recommend adding this reagent in the same order in which the substrate was added.
10. Read the OD of each well at 450 nm within the following 5 min after the addition of the stop solution.

VIII. READING AND RESULT INTERPRETATION

The ODs must be read with an ELISA reader at 450 nm.

If samples and controls are run in duplicate, calculate the arithmetic mean of the OD values obtained.

Validation of the kit:

O.D positive control > 1
O.D negative control < 0,3

Interpretation of the results

The cut off will be calculated using the positive control OD value:

Cut off = OD positive control \times 0,25

Samples should considered:

- ⇒ Negative samples: those with an $OD_{450} \leq$ cut off.
Positive samples: those with an $OD_{450} >$ cut off

Despite the good sensitivity and specificity rates of this test, it is recommended that all positive samples are confirmed using another diagnostic method. The final diagnosis should not be based only on a single test result. In all cases, these results should be corroborated and related to other important information: medical history of the animal, vaccines applied and clinical symptoms.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: