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
TITULO (y) = 54 (e^{4 * 0.681}) = 824

I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the antigen on a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on plate. After washing to eliminate all non-fixed material from the sera sample, we can detect the presence of dog immunoglobulins using a specific peroxidase conjugate. After addition of the substrate a colorimetric reaction will

appear which could be measured by a spectrophotometer.

In this way the presence of colour means the presence of antibodies against the virus in the dog sera, and the absence of colour the absence of specific antibodies.

In our kit is important to remark the use of antigen obtained by the expression of the viral proteins in baculovirus growth in insect cell cultures. This method warranties the total absence of infectivity in the kit.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.
10. The substrate and stop solution, handle with care.

III. STORAGE OF THE REAGENTS:

Stored between +2°C and +8°C.

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use.
- Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

The sera samples needs to be at 1/100 dilution in diluent to be tested (5µl of serum in 0,5 ml of diluent).

VI. PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the kit into 9 parts of distilled or deionized water. When ready this solution remains stable when is stored between +2°C and +8°C.

Preparation of the control sera (+) and (-).

Control sera are ready to use. DO NOT DILUTE.

Preparation of the conjugate: to make immediately before use.

Dilute 1/100 with diluent.

- The necessary and sufficient quantity of conjugate for a complete plate is 110 µl of conjugate in 11 ml of diluent.
- The necessary and sufficient quantity of conjugate for an 8 well strip is 10µl of conjugate in 1 ml of diluent.

Prepare only the quantity needed for each time because the remainder volume has to be rejected.

Shake very well the solution prior to use.

VII. TEST PROCEDURE:

1. Bring all reagents at room temperature before starting the test. (at least. 2 hours)
2. Add 100 µl of positive control serum, to 2 wells of the first row of the plate and 100 µl of the negative control serum to other 2 wells of the first row of the plate. Add 100 µl of sera samples to test, prepared following the previous instructions, on the remainder wells of the plate. (for confirmatory purposes test the samples in duplicate wells)
Seal the plate and **incubate for 10 min at room temperature (18-25°C).**
3. Wash 3 times following the procedure previously described.
4. Add 100 µl of conjugate, prepared following the previous instructions, to each well. Seal the plate and **incubate for 10 min at 18-25°C (room temperature).**
5. Wash 5 times following the described procedure.
6. Add 100 µl of substrate solution, to each well. Keep the plate for **5 min at room temperature.**
7. Add 100 µl of stop solution to each well.
8. Read the absorbances of each well with an ELISA reader at 450 nm within 5 min after the addition of stop solution.

VIII. READING AND RESULT INTERPRETATION

The reading must be done with an ELISA reader at **450 nm**.

Results Validation:

The test is validated when:

- OD 450 Control (+) is higher than 1.2
- OD 450 Control (-) is lower than 0.15

Results Interpretation:

The ratio S/P (sample OD/control+ OD) must be calculated. Samples with S/P lower than 0.15 should be considered negative and samples with S/P ratio higher than 0.15 as should positive to CPV antibodies.

The titre of the sample will be calculated as follow:

$$Y (\text{titre}) = 54 (e^{4x})$$

Where, **e** is the base of natural logarithm (2,718282) and **x** the S/P ratio of the sample.

EXAMPLE:

$$\begin{aligned} \text{OD Control (+)} &\Rightarrow 1,764 \\ \text{OD Control (-)} &\Rightarrow 0,08 \\ \text{OD of sample} &\Rightarrow 1,202 \text{ (S/P = 0.681)} \end{aligned}$$

EI test is valid because (C+) is higher than 1,2 and C (-) is lower than 0,15.

The titre of the sample is:

$$\text{TITRE (y)} = 54e^{4 * 0.681} = 824$$

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in _____ by: