INGEZIM CPV DAS

Prod Ref: 15.CPV.K2

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I. TECHNICAL BASIS

The kit is based on a **Double antibody sandwich enzyme immunoassay**. This technique shows big advantages from the rest of the technique used till now: more sensitivity, more specificity, faster to do and easily automatizable.

Here we describe the technical basis of the assay: and specific monoclonal antibody (Mab) against CPV is coated on a solid support (polystyrene plate). When a sample containing viral antigen is added, the Mab will capture the viral particles. After washing to remove all material not bound to the plate, we add a second specific Mab to CPV labelled with biotin (conjugated I). This Mab will bind

to it, and after incubation with a second conjugate (streptavidin-peroxidase), we could demonstrate the presence of this conjugate (conjugate II) by adding a chromogenic substrate that in presence of the enzyme will produce a colorimetric reaction.

The use of these very specific and stable Mabs, warrants the objectivity, specificity, safety and accuracy of the assay.

The coating Mab as well as the conjugate Mab, are mixtures of very specific Mabs that allow the detection of any of the CPV strains described in Europe.

II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. For each utilisation of the kit, control positive and negative serum must be tested in a systematic way.
9. Use a new tip for each serum sample.

III. STORAGE OF THE REAGENTS:

Store between +2°C and +8°C. (ATTENTION! Stop solution may precipitate at this temperature, if it occurs, warm at 37°C).

IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

1. Throw out the content of the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
 2. Dispense a volume of 300 µl of washing solution on each well.
 3. Shake delicately the plate, avoiding the contamination between wells
 4. Turn over the plate brusquely to empty the wells.
 5. Repeat the process as much times as is indicated on the instructions of the kit.
 6. Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
 7. After the last step of washing shake the plate turned over an absorbent filter paper.
- **ATTENTION!** you are working with biologic material that could contain infectious agents. All rejected material needs to be sterilised. We recommend using a recipient containing NaOH to through out the liquids after washing. After the last step of washing shake the plate turned over an absorbent filter paper.

V. PREPARATION OF SAMPLES:

Homogenise 1 gr of the faeces sample in 2 ml of diluent I. After a good homogenisation of the sample, centrifuge the homogenate at 1000 r.p.m. The supernatant obtained will be used as the assay sample..

VI. PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water(40 ml of concentrate solution and 960 ml of water). Once this solution ready, it remains stable between +2°C and +8°C.

Preparation of the control (+) and (-)

Prepare a 1/100 dilution of the controls provided with the kit in sample diluent (5 µl of control in 500 µl of diluent I).

Preparation of the conjugate I: to make immediately before use.

Dilute 1/100 the needed quantity of conjugate with diluent I (e.g. 10 µl of conjugate I with 1 ml of diluent I is more than enough quantity for a strip of 8 wells). Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

Preparation of the conjugate II: to make immediately before use.

Dilute 1/100 the needed quantity of conjugate II into diluent II (e.g. 10 µl of conjugate II with 1 ml of diluent II, is enough quantity for a strip of 8 wells). Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

Preparation of the substrate solution: to make immediately before use:

Mix 1 volume of the substrate provided in the kit with 9 volumes of substrate buffer (e.i. 100 µl of substrate with 0.9 ml of substrate buffer is the quantity needed for one complete strip of 8 wells). Prepare only the quantity needed for each time because the remainder volume has to be rejected.

VII. TEST PROCEDURE:

1. Take out the plate (or strips to be used) and Bring to room temperature before starting the test.
2. Add 100 µl of positive control, diluted as it's specified in the previous instructions, to one well (e.i. A1) and 100 µl of the negative control (prepared as is specified before) to another well (e.i. B1). Add 100 µl of each sample to test on each remainder wells. We recommend running the samples in duplicate wells. Seal the plate and **incubate for 1h a 37°C**.
3. Throw away the liquid contained in the plate on a recipient containing NaOH 0.1 M. Wash 4 times following the procedure previously described.
4. Add 100 µl of conjugate I, prepared following the previous instructions, to each well. Seal the plate and **incubate for 1 h a t 37°C**.
5. Wash 4 times following the described procedure.
6. Add 100 µl of Conjugate II, prepared following the previous instruction, to each well. **Incubate 15 minutes at room temperature (25°C)**.
7. Wash 4 times more following the described procedure.
8. Add 100 µl of substrate solution, prepared following previous instruction, to each well. Keep the plate for **10 min at room temperature**.
9. Add 100 µl of stop solution to each well.
10. Read the absorbances with a spectrophotometer at 405 nm within 5 min after the addiction of stop solution.

VIII. READING AND RESULT INTERPRETATION

Read the OD of each well with and ELISA reader at **405 mn**.

Kit Validation

$$\frac{\text{OD Positive control}}{\text{OD Negative control}} > 10$$

Results interpretation

CUT OFF = 0,2 x (OD Positive Control)

- ⇒ Samples with OD values, higher than Cut off value, must be considered as positive to CPV antigen.
- ⇒ Samples with OD values, lower than Cut off value, must be considered as negative to CPV antigen.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.
C/ Hnos. García Noblejas, 39
28037- MADRID (SPAIN)
Tlf: +34 91368.05.01/04
Fax: +34 91 408.75.98
E-mail: ingenasa@ingenasa.com
www.ingenasa.es

Distributed in

by: